



Alpha-Leco® section Inorganique

Choisissez votre Gamme

Alpha-Leco XREF

LECO® CS125/225

LECO® CS200/C200/S200/WC200

LECO® CS230

LECO® CS244/344/IR212/232

LECO® CS300

LECO® CS400/C400/S400

LECO® CS444/IR4112/432

LECO® CS600

LECO® DH603

LECO® RH402

LECO® RH404

LECO® RO478

LECO® TC436/TN414/RO416

Vous serez automatiquement redirigé à la section souhaitée

ALPHA-LECO XREF

LECO	ALPHA	DESCRIPTION
145-601	AR-5601	WATER TEMP REGULATOR
156-002	AR-002	"O" RING (UPPER ELECTRODE)
156-004	AR-005	TEFLON BACKUP SEAL
163-012	AR-63012	TUBING INSERT
170-001	AR-0001	SPRING, LOWER ELECTRODE
172-505	AR-2505	TEFLON TUBING, 1/4" OD
175-573	AR-5573	TEAR OFF BLADE
175-574	AR-5574	OPTICAL SWITCH ASSEMBLY
175-590	AR-5590	PUMP REPAIR KIT, SC32/132
175-631	AR-5631	MOTOR BRUSH FOR 774-571
175-632	AR-5632	"O" RING, FLOW CONTROLLER
175-659	AR-5659	PAPER DEFLECTOR V55147
175-673	AR-5673	PRINT HEAD
175-674	AR-5674	TEAR OFF BLADE SPRING
183-001	AR-2174	ALPHASOLVE II, 20 - 30 MESH
188-012	AR-293	SMALL THROAT BRUSH
190-142	AR-142	SCREWS FOR LOWER ELECTRODE, 4 REQUIRED
190-739	AR-739	NYLON SPACER SCREW
194-060	AR-4060	HEX JAM NUT FOR AR-5814
204-001	AR-001	THERMOREGULATOR TC30,TN14,RO16
501-001	AR-1701	COAL STANDARD, 0.4% SULFUR
501-005	AR-1704	COAL STANDARD, 1.0% SULFUR
501-008	AR-008	ALPHACEL II TUNGSTEN,2.2 KG
501-009	AR-1009	SPRING CLIP FOR DUST TRAP
501-010	AR-010	DUST TRAP COMPLETE
501-011	AR-011	DUST TRAP BODY INSERT
501-012	AR-011B	DUST TRAP RESERVOIR
501-013	AR-013	FILTER CLOTHS, PKG. 10
501-014	AR-1705	COAL STANDARD, 1.5% SULFUR
501-020	AR-1708	COAL STANDARD, 3.0% SULFUR
501-022	AR-1709	COAL STANDARD, 3.5% SULFUR
501-024	AR-322	GREY IRON STANDARD 150GM., 3.35%C ,0.041 %S
501-027	AR-204	ASPIRATOR BULB
501-031	AR-031	INDUCTION COIL, HYDROCARBON
501-034	AR-1034	CALCIUM CARBONATE, 12%, 50GM
501-045	AR-1045	ZIRCON CRUCIBLE, CHN 900
501-046	AR-2021	SUCROSE STANDARD, 50G

501-053	AR-1053	ACETANILIDE STANDARD, 10GM
501-057	AR-079G	ALPHANAL, 1 GALLON
501-059	AR-059	TIN CAPSULES, 100 btl
501-060	AR-060	MANGANESE DIOXIDE 30 CC
501-061	AR-1061	ARROWROOT STARCH, 227GM
501-062	AR-061	BOAT PULLER
501-067	AR-067	ROTOMETER BALL ONLY
501-073	AR-073	GRAPHITE POWDER (BTL) (100GMS)
501-076	AR-076	TIN METAL
501-077	AR-077	IRON CHIP
501-078	AR-078	IRON POWDER
501-079	AR-079	ALPHANAL WETTING AGENT, 118 mL
501-081	AR-081	GLASS WOOL
501-082	AR-022	BRUSH
501-083	AR-304	DUCTILE IRON STANDARD, 4.23% C
501-089	AR-507	VALVE ASSEMBLY FOR 516-000
501-095	AR-2814	RESIDUAL OIL, 1.0% S
501-099	AR-2837	LUBE OIL, 4% S
501-101	AR-800	ARROWROOT STARCH SOLUTION
501-105	AR-301	CAST IRON STANDARD, 2.12% C
501-110	AR-110	3/8" I.D. SILICONE TUBING
501-110-019	AR-110-19	SILICONE RUBBER TUBE, 3/8" X 19"
501-112	AR-112	1/4" I.D. NEOPRENE TUBING
501-123	AR-307	TUNGSTEN CARBIDE STD, 100G
501-130	AR-1710	COAL STANDARD, 4.0% SULFUR
501-141	AR-2779	ULTIMATE COAL STANDARD, 50GM
501-144	AR-144	MAGNESIUM OXIDE
501-147	AR-148	COPPER PIN STANDARD, 500PPM
501-148	AR-149	COPPER PIN STANDARD, 3000PPM
501-165	AR-1165	DRIERITE INDICATING 10-20M, 454GM
501-166	AR-166	ASCARITE II 8 - 20 MESH
501-170	AR-170	RARE EARTH COPPER OXIDE
501-171	AR-171	ANHYDRONE, 8-20 MESH
501-179	AR-079P	ALPHANAL, 16 OZ
501-184	AR-1711	COAL STANDARD, 5.0% SULFUR
501-197	AR-197	TRAP MATERIAL 12GM
501-233	AR-233	CERAMIC ADHESIVE ALUMINA, 2OZ
501-234	AR-234	OXYGEN SPECIMEN FILE
501-241	AR-241	SILICONE GREASE 5.3 OZ.

501-263	AR-263	COPPER METAL ACCELERATOR, 3 lb (1.36 kg)
501-264	AR-264	COPPER METAL ACCELERATOR, 0.5 lb (227 g)
501-265	AR-1265	CRUCIBLE HOLDER TRAY
501-280	AR-634	"O" RING
501-291	AR-291	OXYGEN REGULATOR ASSEMBLY
501-369	AR-1727	PROX COAL STD. 2.5% S
501-426	AR-426	COMBUSTION CONTROL, SOLID SAMPLES
501-427	AR-427	COMBUSTION CONTROL, LIQUID SAMPLES
501-428	AR-1721	PROX COAL STD. 0.4% S
501-429	AR-1723	PROX COAL STD. 1.0% S
501-430	AR-1724	PROX COAL STD. 1.5% S
501-431	AR-1726	PROX COAL STD. 2.0% S
501-432	AR-1728	PROX COAL STD. 3.0% S
501-433	AR-1730	PROX COAL STD. 4.0% S
501-438	AR-734	PROX COAL STD. 0.9% S
501-439	AR-000	PARAFIN OIL, 118ML/BTL C, H, N, BTU
501-441	AR-2021	SUCROSE, 50GM
501-457	AR-400	TIN COATED COPPER, 1.3KG
501-458	AR-2822	#2 DIESEL OIL
501-501	AR-870	RING STANDARD, 1# JAR, .031 C
501-502	AR-871	RING STANDARD, 1# JAR, .053% C
501-503	AR-872	RING STANDARD, 1# JAR, .185% C
501-504	AR-873	RING STANDARD, 1# JAR, .468% C
501-505	AR-874	RING STANDARD, 1# JAR, .723% C
501-506	AR-875	RING STANDARD, 1# JAR, .817 %C
501-529	AR-556	1 GM STEEL BALL, 5.7PPM H2
501-539	AR-930	ACTIVATED CHARCOAL, 454GM
501-551	AR-660	N/O STANDARD, 0.5 G, STEEL, 100 BTL
501-552	AR-656	N/O STANDARD, 0.5 G, STEEL, 100 BTL
501-553	AR-669	N/O STANDARD, 0.5 G, STEEL, 100 BTL
501-561	AR-2016	SOY BEAN MEAL SRM
501-563	AR-2025	CORN MEAL SRM
501-571	AR-571	COPPER CAPSULES, 100 BTL
501-587	AR-912	PLATINIZED SILICA
501-592	AR-062	MANGANESE DIOXIDE, 250 CC
501-598	AR-598	NIBBLED NICKEL FLUX, 100GM BTL
501-608	AR-1608	QUARTZ WOOL WADDING, 200 GM
501-609	AR-615	FURNACE REAGENT, 100GM
501-614	AR-1614	SPATULA

501-621	AR-621	COPPER TURNINGS
501-632	AR-632	MOLE SIEVE TRAP 14G
501-636	AR-636	VANADIUM PENTOXIDE
501-640	AR-366	COPPER METAL ACCELERATOR, 25#
501-641	AR-643	COPPER METAL SCRUBBER MATERIAL, 454GM
501-644	AR-644	OXYGEN STANDARDS, 1GM STEEL
501-645	AR-1651	N/O STANDARD, 1GM PINS
501-647	AR-1655	N/O STANDARD, 1.0 G, STEEL, 100 BTL
501-648	AR-2776	ULTIMATE COAL STANDARD, 50GM
501-673	AR-673	IRON CALIBRATION STANDARD, 454GM
501-674	AR-881	PIN STANDARD, .017% C, .003% S
501-675	AR-882	PIN STANDARD, .064% C, .021% S
501-676	AR-890	PIN STANDARD, .099% C, .018% S
501-677	AR-895	PIN STANDARD, .247% C, .021% S
501-678	AR-884	PIN STANDARD, .441% C, .011% S
501-679	AR-885	PIN STANDARD, .805% C, .022% S
501-697	AR-612	PLATINIZED CHARCOAL, 50GM
501-951	AR-1951	NICKEL CALIBRATION WIRE
501-955	AR-1722	PROX COAL STD. .7% S
501-957	AR-1729	PROX COAL STD. 3.5% S
501-959	AR-742	PROX COAL STD. 2.0% S
501-961	AR-304	DUCTILE IRON STANDARD, 4.15% C
501-964	AR-1712	COAL STANDARD, 6.0% SULFUR
501-967	AR-2851	EXTREME PRESSURE OIL
501-968	AR-2852	EXTREME PRESSURE OIL
501-969	AR-2853	EXTREME PRESSURE OIL
502-005	AR-1205	PIPETTE BRUSH
502-007	AR-2189	TUBE CLEANER, 8MM X 4FT, FP428
502-008	AR-570	TIN PLUGS, 150/BTL
502-016	AR-655	N/O STANDARD, 0.5G, STEEL, 100 BTL
502-020	AR-599	ETCHABLE NICKEL FLUX 100GM
502-021	AR-548	HYDROGEN STANDARD, .5PPM H ₂ /100
502-023	AR-473	REAGENT TUBE FUNNEL
502-024	AR-651	0.1GM PINS, 10GM BTL, 29PPM HYDROGEN, TI
502-026	AR-2818	RESIDUAL OIL
502-029		CONTACT ALPHA FOR ENTIRE RANGE
502-030		CONTACT ALPHA FOR ENTIRE RANGE
502-032	AR-078	ACCELERATOR, 454G FE POWDER
502-040	AR-2166	TIN CAPSULE, 6MM OD, 100/BTL

502-042	AR-025	QUARTZ CRUCIBLE, CHN
502-049	AR-049	REAGENT "N" CATALYST
502-060	AR-556	1.0 GM STEEL PIN STD., 5.7PPM HYDROGEN
502-061	AR-555	1.0GM STEEL PIN STD., 2.5PPM HYDROGEN
502-062	AR-4020	SOIL STANDARD, 1% C, 100 G
502-071	AR-2864	HIGH SULFUR RESIDUAL OIL, 118ML 30% S
502-072	AR-655	1/2 GM STEEL PIN, N = 3896 PPM O = 2
502-076	AR-571	COPPER CAPSULES, 100/BLT
502-077	AR-570	TIN PLUGS, 150/BTL
502-085	AR-3085	ZN SULFIDE, 32% S, 50GM
502-091	AR-091	CALCIUM OXALATE STD, 50GM
502-092	AR-2092	EDTA STANDARD, 50GM
502-122	AR-2122	FIRING FUSE WIRE, 2" PKG OF 10
502-134	AR-2134	BENZOIC ACID TABLET 1.45GM, PK 10
502-134	AR-2134V	BENZOIC ACIT TABLET 1.45GM, PK 500
502-154	AR-648	TI PIN STD, 130 PPM HYDROGEN, BTL 100 1/4 G PIINS
502-156	AR-2156	FLUORHIBIT, 454GM
502-167	AR-2167	TIN CAPSULES, 8MM, 100/BTL, FP428
502-173	AR-173	ALPHACEL II - HP, 750G
502-174	AR-2174	ALPHASOLVE II, CO2 ABSORB 20-30, 500GM
502-176	AR-2176	ALPHASOLVE II, CO2 ABSORB 8-20, 500GM
502-177	AR-177	FINE QUARTZ WOOL, 50 GM
502-181	AR-2181	TUNGSTIC OXIDE, GRANULAR, 60GM
502-183	AR-2183	NICKEL CAPSULE, 10.0 mm L X 6.7 mm W, 100/BTL
502-184	AR-502	BENZOIC ACID POWDER, 99%, 10 GM
502-186	AR-2186	TIN FOIL CUPS, 100/BTL, FP428
502-187	AR-2187	TIN CAPSULE, 3.5mm x 9mm, 100/BTL
502-188	AR-2188	ALUMINA OXIDE PELLETS,FP428,200G
502-189	AR-189	COPPER STICKS, 100GM, FP228/428
502-190	AEB-1002	COPPER OXIDE, FINE WIRES, 100GM
502-203	AEB-2031	STEARIC ACID STANDARD, 1GM
502-204	AEB-2034	DINITROBENZOIC ACID STANDARD, 1GM
502-205	AEB-2030	CAFFEINE STANDARD, 1GM
502-206	ASD-2003	SILVER CAPSULES, 3.5 X 5 MM, 100PCS
502-207	AEB-2035	CYSTINE STANDARD, 1GM
502-208	AEB-2036	SULPHANILAMIDE STANDARD, 1GM
502-209	AR-2209	SULPHAMETHAZINE STANDARD, 5GM
502-211	AR-211	GLYCINE 99%, 50 GM
502-213	AR-2213	LEAK DETECTION SOLUTION, 4OZ

502-225	AEB-1036	SILVER TUNGSTATE/MGO, 100GM
502-226	AEB-1009	SILVER VANADATE, 10GM
502-227	ATD-1000	TIN CAPSULES, 3.2mm X 4mm, 500
502-231	AR-673	HI PURITY IRON STD, LOW C & S
502-258	AR-621	COPPER TURNINGS, DEGASSED, 200GM
502-266	AR-205	FIRING FUSE WIRE, 300 DTMS, 0.006" Dia
502-272	AR-2017	CORN GLUTEN, C,H,N,S, 30G
502-273	AR-2018	ALFALFA, C,H,N,S, 30G
502-274	AR-2019	WHEAT FLOUR, C,H,N,S, 30G
502-275	AR-2020	RYE FLOUR, C,H,N,S, 30G
502-276	AR-2026	OAT MEAL, C,H,N,S, 30G
502-277	AR-2027	BARLEY, C,H,N,S. 30G
502-278	AR-2028	RICE FLOUR, C,H,N,S. 30G
502-282	AR-189	COPPER STICKS, 100G AMPULE
502-294	AR-008F	REPLACED BY AR 008 OR AR 173
502-295	AR-621	COPPER TURNINGS, DEGASSED, 200GM
502-297	AR-008B	ALPHACEL II, HP, 13.6KG
502-298	AR-298	SULFAMETHAZINE, 25 GM
502-301	AR-2301	TUNGSTEN OXIDE, 125GM
502-304	AR-2304	COPPER STICKS, 100G PACKED IN ARGON AMPOULE
502-304-500	AR-2304-500	COPPER STICKS, RESEALABLE BTL, 500 GM
502-308	AR-4016	SOIL STANDARD, 2% C & S, 25 G
502-310	AR-2310	STEEL WOOL FILTER MATERIAL, 454GM
502-311	AR-2311	FURNACE REAGENT 1000, 100GM
502-314	AR-177	QUARTZ WOOL, 28.4GM
502-318	AR-4017	SOIL STANDARD, 25 GMS, C=0.5%, S= 0.5%
502-319	AR-4005	ORE STANDARD, 30 GMS, C=1.42%, S=1.41%
502-320	AR-4006	ORE STANDARD, 30 GMS, C=4.40%, S=4.24%
502-321	AR-321	ALPHACAT ACCELERATOR, 250 GM
502-335	AR-2335	NICKEL FOIL CRUCIBLE LINERS, PKG 250
502-338	AR-338	ALPHAGEL CAPSULES, 400 BTL, FOR FP428
502-343	AR-2343	NICKLE BOAT LINER, PKG/100
502-344	AR-2344	NICKEL BASKETS, LOW OXYGEN,1G/100BTL
502-345	AR-2345	NICKEL BASKETS, LOW OXYGEN,1.5G/100BTL
502-351	AR-351	HALOGEN SCRUBBER, 50 G
502-359	AR-2359	FURNACE REAGENT, MGO, 200 GM
502-372	AR-4007	ORE STANDARD, 30 GMS, C=7.62%, S=3.25%
502-374	AR-2374	HEATED PARTICULATE FILTER, PKG OF 50, 25MM
502-382	AR-339	MEDIUM ALPHAGEL CAPS, 400 CAPS PER BTL

502-383	AR-1700	COAL STANDARD, 50 GM, 0.27% S NOM
502-384	AR-1702	COAL STANDARD, 50 GM, 0.74% S NOM
502-385	AR-1704	COAL STANDARD, 50 GM, 0.95% S NOM
502-386	AR-1705	COAL STANDARD, 50 GM, 1.54% S NOM
502-387	AR-1708	COAL STANDARD, 50 GM, 2.03% S NOM
502-388	AR-1711	COAL STANDARD, 50 GM, 2.93% S NOM
502-389	AR-2892	RESIDUAL OIL STANDARD, 0.10%S, 4 OZ (118ML)
502-390	AR-2895	RESIDUAL OIL STANDARD, 0.50%S, 4 OZ (118ML)
502-391	AR-2814	RESIDUAL OIL STANDARD, 1.00%S, 4 OZ (118ML)
502-392	AR-2815	RESIDUAL OIL STANDARD, 2.00%S, 4 OZ (118ML)
502-393	AR-2816	RESIDUAL OIL STANDARD, 3.00%S, 4 OZ (118ML)
502-394	AR-2899	RESIDUAL OIL STANDARD, 4.5%S, 4 OZ (118ML)
502-397	AR-2397	LARGE TIN FOIL CUPS, 100 PER BTL
502-397-400	AR-2397-400	LARGE TIN FOIL CUPS, 400 PER BTL
502-417	AR-2838	SULFUR IN LUBE OIL 0.01%S, 4 OZ (118ML)
502-418	AR-2838	SULFUR IN LUBE OIL, 0.05%S, 4 OZ (118ML)
502-419	AR-2934	SULFUR IN LUBE OIL, 0.10%S, 4 OZ (118ML)
502-420	AR-2833	SULFUR IN LUBE OIL, 0.50%S, 4 OZ (118ML)
502-421	AR-2835	SULFUR IN LUBE OIL, 1.00%S, 4 OZ (118ML)
502-422	AR-2836	SULFUR IN LUBE OIL, 2.00%S, 4 OZ (118ML)
502-423	AR-2837	SULFUR IN LUBE OIL, 4.00%S, 4 OZ (118ML)
502-424	AR-2871	SULFUR IN #2 DIESEL FUEL, 0.01%S, 4 OZ (118ML)
502-425	AR-2873	SULFUR IN #2 DIESEL FUEL, 0.02%S, 4 OZ (118ML)
502-426	AR-2875	SULFUR IN #2 DIESEL FUEL, 0.05%S, 4 OZ (118ML)
502-427	AR-2877	SULFUR IN #2 DIESEL FUEL, 0.10%S, 4 OZ (118ML)
502-428	AR-2823	SULFUR IN #2 DIESEL FUEL, 0.50%S, 4 OZ (118ML)
502-429	AR-2825	SULFUR IN #2 DIESEL FUEL, 1.00%S, 4 OZ (118ML)
502-430	AR-2827	SULFUR IN #2 DIESEL FUEL, 2.00%S, 4 OZ (118ML)
502-431	AR-2829	SULFUR IN #2 DIESEL FUEL, 3.00%S, 4 OZ (118ML)
502-433	AR-1700	COAL STANDARD, 50 GM, 0.29%S NOM
502-434	AR-1702	COAL STANDARD, 50 GM, 0.74%S NOM
502-435	AR-1704	COAL STANDARD, 50 GM, 0.95%S NOM
502-436	AR-1705	COAL STANDARD, 50 GM, 1.54%S NOM
502-437	AR-1708	COAL STANDARD, 50 GM, 2.93%S NOM
502-438	AR-1710	COAL STANDARD, 50 GM, 4.95%S NOM
502-450	AR-450	HEATED PARTICULATE FILTER, PKG OF 50, 3.2MM
502-489	AR-2489	POROUS CRUCIBLE CHN, 10 PKG
502-493	AR-2156	FLUORHIBIT, 454GM
502-810	AR-810	V-CAPS LARGE

502-820	AR-568	NICKEL CAPSULE,12.3mm L X 6.7mm W,100 BTL
502-822	AR-568	NICKEL CAPSULE, 12.3 mm L X 6.7mm W, 100 BTL
502-825	AR-825	TIN CAPSULE 9MM X 21MM, 50 PKG
503-007	AR-007	DETACHABLE FLOAT VALVE
503-020	AR-3020	BORON DISTILLATION APPARATUS
503-031	AR-905	FUNNEL (3 LEG)
503-032	AR-901	SCOOP
503-033	AR-033	SULFUR TRAP, FILLED
503-051	AR-051	OXYGEN BAFFLE FOR 1 1/8" ID TUBE
503-052	AR-052	OXYGEN BAFFLE FOR 1 1/4" ID TUBE
507-010	AR-911H	CATALYST TUBE (HOPCOLITE)
507-011	AR-911E	CATALYST TUBE (UNFILLED)
507-012	AR-012	REDUCER
507-100	AR-7100	CATALYST FURNACE, REBUILT
507-200	AR-7200	CATALYST FURNACE, REBUILT
507-207	AR-911	CATALYST TUBE, FILLED, PT SI
507-207	AR-911R	CATALYST TUBE, FILLED, RARE EARTH
508-300	AR-8300	CATALYST HEATER, REBUILT
508-400	AR-8400	CATALYST HEATER, REBUILT
508-500	AR-8500	CATALYST HEATER, REBUILT
508-600	AR-8600	CATALYST HEATER, REBUILT
511-036	AR-182	4MM I.D. EVACUATED PIN TUBE
511-038	AR-1038	PIN TUBE HOLDER
511-039	AR-181	3MM I.D. EVACUATED PIN TUBE
511-390	AR-183	6MM I.D. EVACUATED PIN TUBE, PKG 10
516-000	AR-906	PURIFYING TRAIN
516-004	AR-504	ROTOMETER
516-005	AR-505	DRYING TOWER
516-006	AR-506	ACID TOWER
516-063	AR-507	NEEDLE VALVE COMPLETE
517-204	AR-204	RUBBER SQUEEZE BULB
518-016	AR-016	BURRETTE
518-024	AR-024	LIGHT DIRECTOR
518-030	AR-017	VESSEL
518-045	AR-045	REWORKED PHOTOCCELL
519-004	AR-004	HYDROCARBON SULFUR TUBE
519-005	AR-900	IGNITOR
519-009	AR-009	PEDESTAL
521-112	AR-114	SILICONE TUBING, 4" X 3/8" I.D.

521-113	AR-113	SILICONE RING SEAL
521-186	AR-187	SOLID STATE RECTIFIER
521-287	AR-287	SUPPRESSOR COIL
521-305	AR-932	DRYING TUBE
528-018	AR-3818	CERAMIC CRUCIBLE
528-028	AR-8028	FILTERING CRUCIBLES, SMALL
528-030	AR-8030	FILTERING CRUCIBLES, LARGE
528-040	AR-040	COVER, 4MM HOLE
528-042	AR-042	POROUS COVERS
528-043	AR-9880	COVERS W/ 10MM HOLE
528-050	AR-3818B	CERAMIC CRUCIBLES, BAGS/500/BOX
528-053	AR-8053	COMBUSTION BOAT, 4" X 1/2"
528-120	AR-8120	LOW SULFUR CRUCIBLE
528-203	AR-8203	COMBUSTION BOAT FOR SC132/32/CR-12, 500 BOX
528-203-250	AR-8203-250	COMBUSTION BOAT FOR SC132/32/CR-12, 250 BOX
528-230	AR-5282	COMBUSTION BOAT, GLAZED FOR LIQUIDS, PKG OF 50
528-313	AR-8313	COMBUSTION TUBE 1 1/8" X 30" OBE
528-314	AR-8314	COMBUSTION TUBE 1 1/4" X 30" OBE
528-331	AR-8331	COMBUSTION TUBE 1 1/8" X 27" TE
528-332	AR-8332	COMBUSTION TUBE 1 1/8" X 30" TE
528-333	AR-8333	COMBUSTION TUBE 1 1/4" X 30" TE
529-007	AR-529	CERAMIC PEDESTAL, 3.05"
529-011	AR-9011	CERAMIC LINER SC32/132/CR-12
529-012	AR-9012	BOAT STOP KIT (2 X AR9012)
529-018	AR-3203	CERAMIC NOZZLE 320-00203 LECOMELT
529-022	AR-9022	BOAT STOP HALF, HONEYCOMB FINE
529-026	AR-026	CRUCIBLE CHN-600, PKG. 50
529-028	AR-9023	BOAT STOP HALF SC132/32
529-031	AR-9031	CONE SUPPORT D
529-042	AR-9042	CERAMIC CRUCIBLE, MAC-400/TGA, SMALL (5)
529-043	AR-9043	CERAMIC COVER, MAC-400/TGA, SMALL (5)
529-044	AR-9044	CONE SUPPORT I
529-046	AR-9046	CERAMIC TUBE TF-10, 1 1/4" X 23" L
529-047	AR-9047	CERAMIC CRUCIBLE, MAC-400 LARGE (5)
529-048	AR-9048	CERAMIC COVER, MAC-400 LARGE (5)
529-054	AR-9023	BOAT STOP HALF SC32/132
529-055	AR-026	CRUCIBLE CHN, BAG OF 50
529-203	AR-8203P	COMBUSTION BOAT SC32/132/CR12
529-204	AR-9204P	COMBUSTION BOAT, S-144DR, BAG OF 50

529-204-250	AR-9204-250	COMBUSTION BOAT, BOX OF 250
532-060	AR-045	REBUILT PHOTOCELL ASSEMBLY
532-065	AR-065	BURRETTE FOR 532-500
532-066	AR-066	VESSEL FOR 532-500
532-500	AR-532	AUTOMATIC SULFUR TITRATOR
533-104	AR-801	STARCH DISPENSING
537-103	AR-103	1/4" O.D. POLYFLOW TUBING
540-140	AR-140	BULB
544-147	AR-4147	LEVELING BOTTLE FOR 572
549-007	AR-907	ARMATURE
549-072	AR-9072	RESTRICTOR
549-076	AR-176	BULB, LIGHT SOURCE FOR 518, 532
550-004	AR-404S	OSCILLATOR TUBE, EIMAC
550-007	AR-5007	45 SEC TIME DELAY TUBE
550-044	AR-044	INSULATOR, PORCELAIN 521/523
550-046	AR-046	RESISTOR, 521/523
550-047	AR-0047	RF CHOKE FOR 521/523
550-048	AR-048	RF CHOKE
550-053	AR-0053	POWER SWITCH FOR 521/523
550-119	AR-019	RIGHT ANGLE BEND
550-120	AR-120	SULFUR COMBUSTION TUBE
550-122	AR-122	CARBON COMBUSTION TUBE
550-123	AR-123	PEDESTAL INSERT
550-124	AR-124	PEDESTAL, 2.50"
550-165	AR-165	SPRING CLIP FOR COMBUSTION TUBE
550-167	AR-550	FURNACE FAN FOR MODELS 521/523, 115VAC
550-182	AR-180	QEGC CRUCIBLE
550-183	AR-179	QUARTZ CRUCIBLE INSERT FOR AR-180
550-184	AR-184	COPPER RINGS
550-187	AR-186	CAPACITOR
563-035	AR-035	HEADER STOPCOCK, 572, 577
563-037	AR-037	REBUILT BURET FOR 577
566-008	AR-006	INDUCTION COIL
568-016	AR-428	"O" RING
572-016	AR-72016	CONTROL VALVE FOR 572
572-029	AR-029	FLOAT VALVE ASSEMBLY
572-031	AR-2031	SPLASHPROOF DEVICE
578-672	AR-672	"O" RING
585-235	AR-200	AL203 WHEEL, 7" X 3/4" H X 1/16"

585-270	AR-270	AL2O3 CUT-OFF WHEEL 7 X 1" X 1/16"
585-359	AR-5359	PIN CUTTER SPACER, FOR 1/2 G PIN
585-360	AR-5360	PIN CUTTER SPACER, FOR 1 G PIN
585-500	AR-2000X	PIN CUTTER, 1 H.P., 230VAC, 50 Hz.
585-600	AR-2000	PIN CUTTER, 1 H.P., 115VAC, 60 Hz.
589-305	AR-9305	1/4" RUBBER SLEEVE
589-364	AR-9364	OVEN FAN, WR12,TC30, DAYTON
589-551	AR-422	"O" RING (AUTO DROP)
589-699	AR-699	LANCE SLEEVE,CS44,IR12/32,AR521
589-859	AR-859	THERMISTORS, WR-12, 750
589-887	AR-887	"O" RING
589-891	AR-9891	TIMER RELAY 750/589
593-100	AR-593	TIMER, SULFUR TITRATOR
598-142	AR-859	THERMISTORS, WR-12, 750
598-157	AR-157	MOISTURE TRAP, FILLED, 750, 598
598-158	AR-158	MOISTURE TRAP, UNFILLED, 750, 598
598-257	AR-257	REPLACED BY 614-591
598-614	AR-613	BULBS - SWITCHES 6.3V (750-100)
601-064	AR-064	QUARTZ HEATER TUBE, CF-10
601-065	AR-1065	PYREX REAGENT TUBE, CF-10
601-158	AR-1158	"O" RING CHN 800
601-159	AR-1159	"O" RING CHN 800
601-132	AR-132	COMBUSTION TUBE CHN 800
601-133	AR-1133	"O" RING (1.0 X 1.25) CHN
601-201	AR-1201	REACTION TUBE DH103
601-390	AR-390	QUICK DISC TUBE FP428
601-397	AR-397	SHIM, COMBUSTION TUBE FP428, PK8
601-437	AR-437	QUARTZ CATALYST TUBE, FP428
601-438	AR-438	COMBUSTION TUBE FOR NP228/428, SN 1 THRU 3445
601-441	AR-441	QUARTZ LANCE TUBE, FP428
601-442	AR-442	"O" RING, FP428, RC612
601-482	AR-482	HONEYCOMB SUPPORT, FP428
601-483	AR-483	PRINTER PAPER, 9 1/2 X 11, 600 SHEETS
601-484	AR-484	PRINTER RIBBON CARTRIDGE
601-490	AR-1490	PELLETIZING FOILER DEVICE
601-504	AR-1504	"O" RING, LOADING BLOCK,FP228
601-505	AR-1505	"O" RING, COMBUSTION TUBE, FP228, RC612
601-506	AR-1506	COMBUSTION TUBE, FP228
601-654	AR-1654	THERMOCOUPLE FOR CHN 2000

601-684	AR-684	"O" RING, LOADING HEAD ASSEMBLY TC600
601-691	AR-1691	.125 CFLEX TUBING, 12"
601-691-030	AR-1691-30	CFLEX TUBE, 30"
601-691-036	AR-1691-36	CFLEX TUBE, 36"
601-862-110	AR-1862	FURNACE, FP-428, 550-650 WATTS
601-920	AR-920	"O" RING
601-939	AR-1939	"O" RING, DOSER, FP428
601-963	AR-1963	TIN CAPSULES, 3.5mm X 4mm, CHN932
601-981	AR-981	SIDE ARM CLEANING BRUSH, FP428
602-615	AR-2615	FURNACE HOUSING W/O END CAPS FP2000
602-802	AR-802	"O" RING, RH 402, RC612
602-803	AR-803	"O" RING, RH 402
602-969	AR-969	5000 DETERMINATIONS PACKAGE FOR FP-428
603-047	AR-63047	URETH TUBING, .061 X .12 OD, RED, PER FT.
603-048	AR-63048	URETH TUBING, .061 X .12 OD, BLUE, PER FT.
603-049	AR-63049	URETH TUBING, .061 X .12 OD, GREEN, PER FT.
603-050	AR-63050	URETH TUBING, .061 X .12 OD, BLACK, PER FT.
603-051	AR-63051	URETH TUBING, .061 X .12 OD, WHITE, PER FT.
603-052	AR-63052	URETH TUBING, .061 X .12 OD, ORANGE, PER FT.
603-350	AR-350	"O" RING EXTRACTION TOOL
603-493	AR-3493	"O" RING
603-706	AR-706	PEDESTAL FOR RH402
603-707	AR-707	REACTION TUBE RH402
603-710	AR-710	QUARTZ CRUCIBLES, RH402, 24PCS
603-910	AR-3910	"O" RING
603-970	AR-3970	"O" RING RH402, UPPER LOAD MECH
603-971	AR-3971	"O" RING RH402, UPPER LOAD MECH
604-112	AR-604	QUARTZ CATALYST TUBE
604-114	AR-923	PEDESTAL ONLY, QUARTZ WITH FITTING
604-163	AR-163	BRASS BRUSH
604-172	AR-172	"O" RING
604-241	AR-4241	REACTION TUBE, HF402/RH402
604-271	AR-271	SS FLUTED BRUSH, DIA SLOTTED
604-373	AR-4373	CAPSULE HOLDER
604-496	AR-4496	TEFLON FLX TUBING .125 X .032W, FT
604-518	AR-518	SS CLEANING BRUSH, UNSLOTTED
604-541	AR-4541	LANCE TUBE CHN 900/932
605-520-087	AR-6055	"O" RING PACKAGE,CS200/300/400
605-814	AR-5814	PEDESTAL ADAPTER ASSEMBLY

605-815	AR-5815	PEDESTAL ADAPTER, CS144 - CS444
605-816	AR-816	CERAMIC PEDESTAL, CS144-444 SiN
605-817	AR-817	PEDESTAL ADAPTER, CS125
605-818	AR-818	PEDESTAL ADAPTER 606/816, EC12,HF10/20
605-828	AR-828	"O" RING, FP428, SN 1 THRU 3445
605-830	AR-830	"O" RING, CHN
605-875	AR-5875	FURNACE RED/OX TUBE, CHN900
606-026	AR-6026	GRAPHITE CRUCIBLE, PKG 10, MICRO OXYGEN
606-202	AR-202	FURNACE FAN FOR FP2000 SERIES
606-224	AR-6224	THERMOCOUPLE, SC432 SERIES
606-286	AR-563	SHORT STRAP, SC432
606-308	AR-6308	STEPPED CERAMIC BOAT STOP, SC432
606-309	AR-6309	OUTER COMBUSTION TUBE, SC432
606-310	AR-1477	"O" RING COMB TUBE, SC432
606-312	AR-1476	"O" RING COMB TUBE, SC432
606-313	AR-6313	INNER COMBUSTION TUBE, SC432
606-316	AR-6316	LANCE TUBE ASSEMBLY, SC432
606-317	AR-6317	STAINLESS PURGE TUBE
606-318	AR-6318	HEATING ELEMENTS (3), SC432
606-319	AR-6319	LEAD WIRE BRAIDED 6.5, SC-432
606-320	AR-6320	BRAID CLAMP FOR SC432 ELEMENT
606-327	AR-6327	GAS DRYING TUBE, SC-432
606-333	AR-6333	"O" RING, LANCE, SC-432
606-340	AR-6340	DRILLED LANCE TUBE, FP-2000
606-403	AR-403	COLLECTION TUBE BRUSH, SC432
606-444	AR-6644	#2 STEPPED CERAMIC BOAT STOP
606-547	AR-6547	BOAT PULLER, SC432 AUTOLOADER
606-550	AR-6550	FINGER ARM FOR AR-6547
606-601	AR-6601	HEATING ELEMENTS (4) CS432,CS144
606-601-100	AR-6601-100	HEATING ELEMENT KIT
606-602	AR-6602	FURNACE END CAP, FP-2000
606-316	AR-6316	LANCE TUBE FOR SC432
606-655	AR-6655	GLASS PURIFIER TUBE, SC432
606-747	AR-6747	2 HOLE LANCE TUBE, FP-2000, 245MM
608-065	AR-8065	LANCE TUBE TRUSPEC MICRO OXYGEN
608-339	AR-8339	COMB TUBE, OLD STYLE, CHN1000 14MM OD X 200MM
608-342	AR-342	COMBUSTION TUBE, QUARTZ, CHN1000
608-343	AR-8343	FURNACE ELEMENT, CHN1000
608-344	AR-344	COMBUSTION TUBE, OUTER, CHN1000

608-378	AR-1476	"O" RING, CHN 1000
608-379	AR-8379	QUARTZ WOOL STRIPS, 15", PKG 10
608-505	AR-8505	COMBUSTION TUBE, CERAMIC, CHN1000
608-866	AR-866	CFLEX TUBING .187" PER FOOT
609-114	AR-9114	"O" RING, SIDE ARM, SN 3446 UP
609-116	AR-9116	COMBUSTION TUBE, FP428, SN 3446 UP
609-129	AR-9129	LANCE, STRAIGHT
609-170	AR-9170	QUICK CAP KIT
609-171	AR-9171	PLUNGER FOR CELLULOSIC TUBES
609-172	AR-9172	TUBING CELLULOSIC, .210ID, 20/PK
609-271	AR-9271	RADIAL SHAFT SEAL, FP2000 AUTOLOADER
609-394	AR-394	HEAT SHIELD FOR FP/NS2000
609-474	AR-9474	FINGER ARM FOR AR-1207 SHAFT
609-728	AR-9728	SHAFT SLEEVE BEARING, FP2000 AUTOLOADER
610-735	AR-735	FURNACE BRUSH, CS-444, DIA SLOTTED
611-207	AR-1207	PULLER SHAFT ASSEMBLY WITH FINGER
611-255	AR-9171	PLUNGER ASSEMBLY
611-351-181	AR691	ELECTRODE INSERT
611-351-182	AR-691	ELECTRODE INSERT
611-351-223	AR-1351	COMBUSTION TUBE, HD
611-431	AR-1431	"O" RING, 132
611-476	AR-1476	"O" RING, INNER COMBUSTION TUBE
611-477	AR-1477	"O" RING, OUTER COMBUSTION TUBE
611-790	AR-790	"O" RING
612-340	AR2340	"O" RING, THERMISTOR BLOCK, TRUSPEC
613-563	AR-563	SHORT STRAP, SC432 PKG 5
614-591	AR-257	CABINET FAN, 110V
614-822-101	AR-4101	SMALL BOAT LOOP
614-822-102	AR-4181	NICKEL BOAT, SMALL FOR AMA 254 Hg
614-822-113	AR-4113	LARGE BOAT LOOP
614-822-114	AR-4114	NICKEL BOAT, LARGE FOR AMA 254 Hg
614-961	AR-4961	CRUCIBLE, FP528, TRUSPEC, PKG 10
614-962	AR-4962	LANCE, FP528
616-063	AR-6063	BOAT PULLER, SULFUR/CARBON
616-068	AR-6068	"O" RING, INNER COMBUSTION TUBE S/SC144
616-070	AR-6070	INNER COMBUSTION TUBE S/SC-144
616-086	AR-6086	STEPPED BOAT STOP, S/SC-144
616-089	AR-6089	OUTER COMBUSTION TUBE, S/SC-144
616-108	AR-6108	COMBUSTION TUBE, FP528

616-138	AR-6138	"O" RING, FP-528
616-139	AR-6139	"O" RING, FP-528
616-146	AR-6146	LANCE, FP528
616-504	AR-6504	PARTICLE FILTER, FP528
616-622	AR-6622	MICRO SWITCH CN 2000 SWITCH ONLY, NO LEADS
616-632	AR-6632	LANCE, SC-432 & FP-2000, 253 mm
616-635	AR-394	HEAT SHIELD FOR FP/NS2000, REPLACES 609-394
617-083	AR-7083	"O" RING, IN SCRUBBER CAT HEATER
617-136	AR-7136	"O" RING, LOADING HEAD, TRUSPEC
617-440	AR-440	"O" RING, FILTER ASSMBY TC600
617-441	AR-7441	MICOR SCREEN FILTER, TCH600
617-486	AR-7486	FITTING REMOVAL TOOL
617-548	AR-7548	FILTER TUBE, TC500, RC612
617-605	AR-7605	LARGE CRUCIBLE FOR LIQUIDS, PKG 10
617-716	AR-7716	QUARTZ FILTER TUBE, TC500
618-376	AR-376	ELECTRODE INSERT
618-482	AR-691	ELECTRODE INSERT
618-603	AR-603	CATALYST TUBE, 13.0"
618-657	AR-6186	FTG ADAPTER, RC612
618-809	AR-8809	REAGENT TUBE, TC600
619-065	AR-9065	U-SHAPED COMBUSTION TUBE FOR TRUSPEC
619-154	AR-9154	CATALYST HEATER TUBE FOR TRUSPEC
619-192	AR-9192	"O" RING, BOTTOM LOAND HEAD, TRUSPEC
619-233	AR-9233	THERMOCOUPLE FOR TRUSPEC TYPE K
619-268	AR-9268	HIGH SENSITIVITY REAGENT TUBE FOR TRUSPEC
619-316	AR-9316	LARGE VOLUME REAGENT TUBE FOR TRUSPEC
619-425	AR-9425	STRAIGHT REAGENT TUBE, TRUSPEC
619-455	AR-9455	LARGE REAGENT TUBE
619-566	AR-9566	BRAIDED FURNACE STRAP, PKG 5
619-590-861	AR-861	"O" RING CS744/844
619-591-385	AR-1385	"O" RING, COMBUSTION TUBE CS844
619-591-683	AR-683	COMBUSTION TUBE - (requires AR683k Kit)
619-591-699	AR-1699	PARTICLE FILTER
619-774	AR-9774	LANCE, TRUSPEC CHNS MICRO
619-776	AR-9776	REAGENT TUBE CHNS TRUSPEC
619-777	AR-9777	DIRECTING MIDDLE TUBE, CHN TRUSPEC
619-803	AR-3302	UPPER ELECTRODE BRUSH, 21 MM DIA
619-880	AR-9880	10 mm HOLE COVER, BOX OF 1000
619-880-100	AR-09880-100	10 mm HOLE COVER, PKG OF 100

619-895	AR-619	GRAPHITE CRUCIBLE, RHN600, PKG OF 1000
619-895-100	AR-619-100	GRAPHITE CRUCIBLE, RHN600, PKG OF 100
620-985	AR3930	UPPER ELECTRODE REPAIR
620-814	AR-691	ELECTRODE INSERT
621-331	AR-9047	TGA CRUCIBLE, PKG OF 5
622-001-026	AR1026	GLASS TUBE 628 SERIES ANALYZERS
625-067	AR5067	HEATING ELEMENTS, PKG 6
625-068	AR-25068	ELEMENT CLIPS, PKG 12
625-401-317	AR1317	PURGE TUBE, TRUSPEC, S.S. 97 MM, 3 HOLES
625-401-352	AR1352	LANCE TUBE, TRUSPEC 225 MM
625-501	AR-6255	THERMOCOUPLE, RC 612
625-501-358	AR1358	THERMOCOUPLE, SC632
625-505-125	AR-5125	CATALYST HEATER TUBE RC612
625-505-138	AR-5138	HEATING ELEMENT ASSMBY RC612
625-505-271	AR-5271	KNURLED ADAPTER, RC612
625-505-168	AR-5168	COMBUSTION TUBE, RC612
625-510-402	AR-5104	COMBUSTION TUBE, RC612
625-601-554	AR-1554	COMBUSTION TUBE, TRUMAC
625-601-555	AR-1555	INNER COMBUSTION TUBE, TRUMAC
625-511-154	AR-511154	LOWER CRUCIBLE CHN, PKG 5
625-602-187	AR-602	CERAMIC LANCE TUBE, 4 HOLES TRUMAC 291.5MM
625-710-392	AR0392	REAGENT TUBE FOR OHN836
625-711-116	AR1116	REAGENT TUBE OHN836
625-711-117	AR1117	REAGENCY TUBE OHN836
660-010-130	AR-9233	THERMOCOUPLE FOR TRUSPEC TYPE K
701-030	APT-1030	TORCH ASSEMBLY, ICP-2000 LECO
702-557	AR-2557	"O" RING, TRUSPEC
750-301	AR-933	DRYING TUVE, 9.0"
751-300	AR-1300	ELECTRONIC BALANCE, 4TH PLACE
760-034	AR-034	GRAPHITE CRUCIBLES
760-072	AR-072	28 VDC BULBS
760-138	AR-138	TWEEZERS, STRAIGHT
760-139	AR-766	TWEEZERS, CURVED DELUX
760-144	AR-022B	BRUSH, 90" BEND
760-145	AR-0145	POWDER FUNNEL FOR GRAPHITE CRUCIBLES
760-193	AR-193	METAL SCREEN
760-224	AR-224	"O" RING
760-225	AR-425	"O" RING (LOADING HEAD PLUNGER)
760-226	AR-428	"O" RING

760-293	AR-293	BRUSH, CLEANING, SMALL
760-304	AR-0304	TEFLON WASHER
760-316	AR-316	WASHER
760-317	AR-317	GLASS TUBE, 6.64"
760-333	AR-333	DRYING TUBE, 8"
760-335	AR-424	"O" RING
760-359	AR-359	"O" RING (QUICK DISCONNECT)
760-414	AR-414	GRAPHITE CRUCIBLES
760-841	AR-194	"O" RING
761-047	AR-047	3" DRYING TUBE
761-133	AR-133	DRYING TUBE, 17.125"
761-190	AR-190	DRYING TUBE, 11.75"
761-375	AR-3931	ELECTRODE REPAIR SERVICE
761-376	AR-3931	ELECTRODE REPAIR SERVICE
761-383	AR-293H	SMALL BRUSH HANDLE ONLY
761-492	AR-492	"O" RING, FP428
761-618	AR-618	TRAP REPAIR CHARGE/WR-12
761-651	AR-1609	"O" RING, LANCE, FP/CHN
761-739	AR-639	TIN PELLETS
761-747	AR-611	SCHUTZE'S REAGENT, 100GM
761-918	AR-918	"O" RING
761-929	AR-929	CRUCIBLE TONGS FOR AR-3818
762-023	AR-023	NUT ADAPTOR SEAL FOR CS44/IR12
762-056	AR-056	COMBUSTION TUBE
762-057	AR-416	"O" RING (DUST FILTER)
762-058	AR-058	"O" RING
762-059	AR-614	"O" RING (CUP ASSEMBLY)
762-060	AR-417	"O" RING BACKUP
762-062	AR-913	OXYGEN LANCE
762-066	AR-7620	H V RECTIFIER,8 kV 1 AMP
762-192	AR-192	IR SO2 CELL, CS44/IR32/33
762-357	AR-357	"O" RING (UPPER ADAPTER)
762-501	AR-2501	HANDLE FOR REBUILT 521
762-523	AR-523	OVEN FAN REPLACES HOWARD FAN
762-621	AR-2621	IR CO2 CELL, CS44/IR12
762-627	AR-2627	BISCUIT FAN, CS44
762-695	AR-695	TIN METAL
762-697	AR-697	LARGE THROAT BRUSH
762-741	AR 650	HYDROGEN STANDARD, 1/4 G X 25G

762-747	AR-546	HYDROGEN STANDARD, 1 .1 PPM H2/100
762-823	AR-697H	BRUSH WITH HANDLE
762-936	AR-936	GLASS SCOOP FOR V2O5, EXTRA LARGE
763-026	AR-027B	ALPHACEL III TUNGSTEN 30# - 20 + 40
763-027	AR-266B	ALPHACEL TUNGSTEN 30# - 12 + 20
763-028-500	AR-027V	ALPHACEL III, TUNGSTEN, 227.27KG
763-029	AR-509	NICKEL BASKETS, 1.5 GM
763-065	AR-508	NICKEL BASKETS, 1.0 GM
763-191	AR-103	TUBING POLYFLO, 1/4" O.D.
763-213	AR-213	GRAPHITE CRUCIBLE FOR RH-1 E
763-261	AR-261	FILTER PAPER, GLASS TUBE
763-263	AR-266	ALPHACEL 5# - 12 + 20
763-265	AR-265	GLASS WOOL, 227GM .
763-266	AR-027	ALPHACEL TUNGSTEN, - 20 + 40
763-295	AR-295	THERMOSTAT, UPPER ELECTRODE
763-379	AR-3379	"O" RING, FP428
763-467	AR-467	IRON CHIP, 25#
763-873	AR-3873	ISOLATION BOX BULB, 28VDC
763-930	AR-930	ACTIVATED CHARCOAL, 454G
764-015	AR-015	BULB FOR 532-500, 28VDC
764-078	AR-4078	1 GRAM CALIBATION WEIGHT
764-098	AR-098	CELLULOSE MATERIAL, 454GM
764-216	AR-4216	HELIUM REGULATOR ASSEMBLY
764-242	AR-242	TIN FLUX, 454GM BTL
764-330	AR-331S	GRAPHITE CRUCIBLE, RH1-EN
764-490	AR-490	PEDESTAL, 2-3/4"
764-502	AR-4544	ROTRON FURNACE FAN
764-545	AR-2832	LUBE OIL, 4 OZ, 0.25% SULFUR
764-546	AR-2835	LUBE OIL, 4 OZ., 1.00% SULFUR
764-547	AR-2836	LUBE OIL, 4 OZ., 2.00% SULFUR
764-640	AR-334	100 PF CAPACITOR
764-641	AR-335	200 PF CAPACITOR
764-720	AR-064	QUARTZ DISCONNECT TUBE, CF-10
764-775	AR-775	WR-12 TRAP MATERIAL
765-068	AR-5068	IR OXYGEN CELL FOR RO17
765-523	AR-523	OVEN FAN, HOWARD KIT
765-624	AR-624	TEFLON TUBING, 1/4" O.D./FT.
765-717	AR-917	OXYGEN LANCE (WR 12)
765-746	AR-746	ROTOMETER, WR12

765-818	AR-5818	OPTIC ISOLATION CABLE
765-933	AR-5933	LUBRICOOL COOLING AGENT
765-976	AR-5976	"O" RING, AC AC100, FP428
766-036	AR-6036	PNEUMATIC (AIR) REGULATOR ASSEMBLY
766-053	AR-766	DELUXE CRUCIBLE TWEEZERS, CURVED
766-363	AR-363	GRAPHITE CRUCIBLES
766-557	AR-6557	"O" RING
766-646	AR-6646	AIR FILTER COATING
766-823	AR-823	MALE CLIPLOCK FITTING
767-243	AR-243	FEMALE CLIPLOCK FITTING
767-277	AR-277	GRAPHITE CRUCIBLES
767-385	AR-085	QUARTZ WOOL, COARSE, 28.4GM
767-500	AR-914	LANCE, IR33
767-518	AR-7518	RETAINING WIRE CLIP SPRING
767-541	AR-541	DRYING TUBE
767-542	AR-542	SCOOP FOR TUNGSTEN, 0.01cc, SMALL
767-543	AR-3931	ELECTRODE REPAIR SERVICE
767-544	AR-3931	ELECTRODE REPAIR SERVICE
767-545	AR-3931	ELECTRODE REPAIR SERVICE
767-546	AR-3931	ELECTRODE REPAIR SERVICE
767-547	AR-3931	ELECTRODE REPAIR SERVICE
767-548	AR-3931	ELECTRODE REPAIR SERVICE
767-549	AR-3931	ELECTRODE REPAIR SERVICE
767-757	AR-1706	COAL STANDARD, 2.0% SULFUR
767-758	AR-1707	COAL STANDARD, 2.5% SULFUR
767-759	AR-720	COKE STANDARD, 0.85% SULFUR
767-838	AR-838	"O" RING
768-774	AR-043S	COVERS, PKG.OF 50-SAME AS AR-043
768-798	AR-8798	CLEAN OUT PLUG FOR SC32/132/CR
768-944	AR-944	"O" RING
768-980	AR-980	PARTICLE FILTER - SC-32/132/RC612
769-335	AR-469	SS RELAY, HF10, 44F 6848 240 D45
769-473	AR-3931	ELECTRODE REPAIR SERVICE
769-520	AR-9520	GRAPHITE CRUCIBLES/RH-2
769-608	AR-608	HALOGEN TRAP MATERIAL, 200GM
769-610	AR-610	HALOGEN TRAP MATERIAL, 225GM
769-743	AR-3931	ELECTRODE REPAIR SERVICE
769-761	AR-761	CRUCIBLES, RH-3, PKG. 100
769-802	AR-3931	ELECTRODE REPAIR SERVICE

769-821	AR-9821	MOLECULAR SIEVE FOR RH-2
769-897	AR-897	GRAPHITE PEDESTAL RH-2 & RH-3
769-898	AR-898	LOWER PEDESTAL
769-994	AR-994	GRAPHITE CRUCIBLE FOR RH-1
770-247	AR-247	FILTER FOR SC32/132, PKG 4
770-436	AR-523	OVEN FAN WITH ADAPTER
770-718	AR-0718	PRINTER
770-764	AR-764	PRINTER PAPER
770-778	AR-778	RIBBON SPOOL ASSEMBLY
770-779	AR-7707	PRINT HEAD WITH CABLE
771-015	AR-1015	MANUAL WT COMP CONTROL CS46
771-030	AR-030	BUBBLE TUBE
771-040	AR-3913	"O" RING
771-275	AR-1275	PNEUMATIC CYLINDER, CS46/IR112
771-354	AR-354	OSCILLATOR TUBE CS-144-444
771-367	AR1367	QUARTZ CRUCIBLES, DC12, PKG 50
771-373	AR-373	DRYING TUBE
771-401	AR-401	REACTION TUBE FOR RH-2
771-407	AR-407	"O" RING
771-425	AR-1425	POWER TRIODE RH-2, 7.5V
771-446	AR-429	"O" RING
771-606	AR-606	PEDESTAL CS-46 & EC-12
771-883	AR-083	BUBBLE TUBE (EC 12)
771-991	AR-991	"O" RING O28
772-005	AR-760	CHART PAPER ROLL, RH2/3
772-269	AR-269	"O" RING
772-381	AR-2381	LANCE CLEANING WIRE
772-402	AR-2402	LANCE SC-32/132/CR-12
772-403	AR-6403	DRILLED LANCE TUBE, SC432
772-446	AR-446	GLASS VIEWING SHIELD, SC32/132
772-449	AR-449	ROD, EXTENSION, GLASS SHIELD
772-456	AR-456	FLOW CONTROLLER ASSY, REPAIR
772-457	AR-456	FLOW CONTROLLER W/O CARD, REPAIR
772-460	AR-2460	PRESSURE SENSOR
772-469	AR-469	REPLACED BY 781-165 (AR-469)
772-508	AR-962	ELEMENT FOR SC32/132/CR12
772-520	AR-520	"O" RING SC32, CS244, 344, RC612
772-522	AR-7073	THERMOCOUPLE ONLY
772-560	AR-2560	CHOPPER MOTOR

772-568	AR-233Q	CERAMIC ADHESIVE ALUMINA, 32OZ
772-738	AR-738	"O" RING
772-739	AR-544	CATALYST TUBE (10.75")
772-747	AR-2747	FAN FOR CS244
772-795	AR-795	OVEN FAN, CS244 & TC136
772-905	AR-057	COMBUSTION TUBE, CS 144-444
772-906	AR-2906	COMBUSTION TUBE RING HOLDER
772-910	AR-2910	"O" RING, PISTON SEAL, FP228
773-066	AR-3066	STAND OFF INSULATOR
773-076	AR-3076	INDUCTION COIL, CS46/HF-768
773-106	AR-106	DELIVERY TUBE
773-136	AR-3136	FITTING, SS FOR FP428
773-140	AR-3140	INNER METAL FILTER, CS-46
773-141	AR-320	"O" RING
773-311	AR-311	SST FILTER
773-314	AR-3314	UNFLUTED CLEANING BRUSH
773-316	AR-789	LANCE TUBE,CS125
773-320	AR-320	"O" RING
773-561	AR-561	SS CLEANING BRUSH, STRAIGHT SLOTTED
773-579	AR-579	METAL SCOOP
773-662	AR-3662	POLYFLO FLEX TUBING, CHN
773-913	AR-3913	"O" RING
773-928	AR-3928	DELIVERY TUBE FOR SC32/132
773-931	AR-9931	TEFLON DISCONNECT BLOCK, SC132
773-932	AR-3932	COMBUSTION TUBE REPAIR SERVICE
774-031	AR-4031	QUICK DISCON TUBE 1.38" OD
774-174	AR-4174	PURGE TUBE
774-204	AR-4204	AC CRUCIBLES, SST, PKG OF 10
774-205	AR-205	FUSE WIRE, 300 DETERMINATIONS, 0.006" DIA
774-208	AR-208	BENZOIC ACID TABLETS, VIAL OF 10
774-211	AR-4211	"O" RING, 2.6 X 2.875
774-253	AR-4253	"O" RING
774-283	AR-4283	"O" RING
774-425	AR-432	"O" RING, CLOSURE BLOCK
774-460	AR-460	THERMOCOUPLE ASSY SC32/132/CR12
774-489	AR-489	"O" RING, UPPER ELECTRODE
774-570	AR-4570	VACUUM PUMP FOR EARLY SC32 110VAC
774-571	AR-4571	VACUUM PUMP FOR SC32/132,24VDC
774-613	AR-4613	"O" RING, HF400 PINCH VALVE

774-638	AR-638	DISCO TUBE, 6" STRAIGHT, RO116
774-649	AR-653	CLEANING BRUSH WITH HANDLE
774-651	AR-652	BRUSH HANDLE ONLY
774-652	AR-652	BRUSH FOR TC-136
774-690	AR-690	KNURLED RING COLLAR,EF100/EF500
774-691	AR-691	REPLACED BY 618-482
774-692	AR-692	LOWER ELECTRODE EXCHANGE/REPAIR ONLY
774-693	AR-693	LOWER ELECTRODE PLATE, REPAIR ONLY AUTOLOAD
774-809	AR-809	INDUCTION COIL HF10, S/OIL
774-841	AR-841	THERMISTOR SENSOR, LOWER DRAIN
775-288	AR-3403	BENZOIC ACID POWDER, 454 GM
775-306	AR-5306	SECONDARY FILTER
775-312	AR-312	"O" RING, LOWER ELECTRODE
775-313	AR-313	TEFLON BACKUP RING
775-314	AR-314	"O" RING, PISTON BLOCK
775-330	AR-5330	FOUR WAY VALVE
775-355	AR-355	CATALYST TUBE, 10.75"
775-416	AR-293	CLEANING BRUSH, TC-136, TN-14
775-416	AR-293H	BRUSH HANDLE
775-431	AR-431	INNER GRAPHITE CRUCIBLE 1000/PK
775-433	AR-433	OUTER GRAPHITE CRUCIBLE 50/PK
775-446	AR-5446	NEW UPPER ELECTRODE
775-446	AR-3930	ELECTRODE INSERT REPAIR SERVICE
775-601	AR-601	BUBBLE TUBE, LEFT
775-625	AR-625	MOMENTARY PUSH SWITCH
775-677	AR-677	"O" RING
775-722	AR-456	FLOW CONTROLLER, REPAIR SERVICE
775-814	AR-402	REACTION TUBE, THICK WALL, RH2/3
775-859	AR-5859	WIRE CLEANING BRUSH, SC32/132
775-880	AR-880	QUICK DISCO TUBE, 4", STRAIGHT, FP228
775-892	AR-431P	INNER GRAPHITE CRUCIBLE, PKG 100
776-017	AR-6017	"O" RING, BALLAST PLATE, TRUSPEC
776-095	AR-095	SILICONE TUBING DUST COLL.
776-129	AR-129	INDUCTION COIL, CS-144-444
776-234	AR-6234	10 MICRON FILTER, 1GPM
776-247	AR-6247	GRAPHITE CRUCIBLES
776-265	AR-6265	ADJUSTABLE PEDESTAL ASSBLY, 923/4060/6922
776-279	AR-279	CATALYST HEATER TUBE
776-283	AR-283	QUARTZ LANCE CHN-600/800

776-284	AR-284	"O" RING, COMB TUBE FP228
776-286	AR-286	COMBUSTION TUBE CHN 600
776-685	AR-6685	TYNA-MITE SOLENOID
776-922	AR-6922	SHIELD & COLLAR FOR PEDESTAL
776-923	AR-923	QUARTZ PEDESTAL POST CS144-444
776-931	AR-606	PEDESTAL HF10/20
776-962	AR-962	ELEMENT FOR SC32/132/CR12
776-978	AR-978	GLASS SCOOP, MAC 400
777-060	AR-7060	QUARTZ COMBUSTION TUBE, SC32/132
777-064	AR-928	PYREX DELIVERY TUBE FOR AR-7060
777-066	AR-7066	QUARTZ LINER FOR AR-7060
777-073	AR-7073	THERMOCOUPLE
778-114	AR-8114	"O" RING 111
778-116	AR-8116	"O" RING 256, BALLAST, FP228
778-120	AR-781	NEOPRENE GASKET, BALLAST, FP228
778-379	AR-8380	COMBUSTION TUBE INSULATOR, FP528
778-492	AR-456	FLOW CONTROLLER, REPAIR SERVICE
778-596	AR-596	SST FILTER
778-674	AR-674	GRAPHITE CRUCIBLE, 2GM
778-891	AR-8891	CRUCIBLE, QUARTZ, LARGE, EA
778-892	AR-8892	COVER, QUARTZ, LARGE, EA
780-044	AR-609	ASSY SAMPLE PULLER
780-047	AR-1202	REACTION TUBE FOR DH102
780-093	AR-093	QUARTZ CRUCIBLE, MAC 400
780-094	AR-094	QUARTZ CRUCIBLE COVER, MAC 400
780-379	AR-379	QUARTZ WINDOW
780-831	AR-831	"O" RING, CHN 1000 LANCE
780-853	AR-853	THERMISTOR, LEVEL
780-890	AR-7808	LARGE GRAPHITE CRUCIBLE, 1000/pkg
780-892	AR-0892	GRAPHITE COVERS FOR AR-7808,1000/pkg
780-899	AR-899	HONEYCOMB SUPPORT, 1.25W X 1.5L FP
780-911	AR-910	CATALYST HEATER TUBE
780-927	AR-927	PRINTER PAPER ROLL
780-928	AR-926	PRINTER RIBBON
780-953	AR-922	PEDESTAL CS-125
781-098	AR-1098	THERMISTOR, HIGH WATER
781-101	AR-101	THERMISTOR, LOW WATER
781-165	AR-469	SS RELAY TC136 44F 6848 240 D45
781-273	AR-028	CERAMIC CRUCIBLE TRAY

781-335	AR-1335	QUARTZ BOAT FOR RC-412, RC612
781-371	AR-371	FAN FOR CS244, 240V, 50/60Hz
781-424	AR-1424	"O" RING
781-585	AR-1585	OVER TEMP THERMOCOUPLE, RC612
781-671	AR-671	INK JET PRINTER PAPER
781-698	AR-1698	CATALYST TUBE RH402,QUARTZ,10.75"
781-883	AR-1883	INK JET PAPER 500
781-902	AR-902	INK JET CARTRIDGE
782-004	AR-928	DELIVERY TUBE SC132/32/CR-12
782-010	AR-7060	QUARTZ COMBUSTION TUBE
782-038	AR-2038	COMBUSTION TUBE FOR RC-412
782-059	AR-2059	NICKEL BOAT LINERS, PKG/10 FOR AR1335
782-079-110	AR-9023	CERAMIC BOAT STOP, HI FLOW
782-079-115	AR-9110	CERAMIC BOAT STOP, MEDIUM FLOW
782-079-120	AR-9022	CERAMIC BOAT STOP, LOW FLOW
782-338	AR-2338	MOISTURE REMOVAL TUBE, INDICATING
782-428	AR-2428	QUARTZ CATALYST TUBE, RC-412
782-534	AR-3930	ELECTRODE INSERT REPAIR SERVICE
782-601	AR-2601	VACUUM CLEANER
782-703	AR-7203	HI PURITY OUTER GRAPHITE CRUCIBLE 1000/pkg
782-703-100	AR-2703P	HI PURITY OUTER GRAPHITE CRUCIBLE 100/pkg
782-719	AR-782P	HI TEMP GRAPHITE CRUCIBLE, 200/pkg
782-720	AR-7827	HI TEMP GRAPHITE CRUCIBLE, BOX OF 1000
782-720S	AR-720S	GRAPHITE HI-DENSITY HI-TEMP CRUCIBLE
782-720SHD	AR-782HD	GRAPHITE HI-DENSITY HI-TEMP CRUCIBLE
782-721	AR-721	HI-TEMP ELECTRODE INSERT
782-795HD	AR-2795	INNER CRUCIBLE FOR AR-2703, 1000/pkg
782-835	AR-3930	ELECTRODE INSERT REPAIR SERVICE
782-877	AR-877	TRIM STEP BRUSH FOR 782-918
782-880	AR-2880	VACUUM FILTER FOR CS444
782-918	AR-877H	HANDLE ONLY FOR AR-877
782-938	AR-938	COMPOSITE BRUSH FOR 782-942
782-942	AR-938H	HANDLE ONLY FOR AR-938
783-177	AR-962	HEATING ELEMENT FOR SC32/132, CR12
783-189	AR-3189	"O" RING, CHN1000
783-190	AR-3897	"O" RING, CHN1000 REPLACED BY 783-897
783-193	AR-8343	CHN-1000 FURNACE ELEMENT WITH INSULATION
783-274	AR-278	LANCE ONLY WITHOUT METAL RING
783-278	AR-278	FURNACE LANCE, CHN1000

783-299	AR-2038	COMBUSTION TUBE FOR RC-412
783-301	AR-3301	BRUSH END, .38DIA.
783-302	AR-3302	BRUSH END MODIFIED
783-515	AR-515	"O" RING
783-897	AR-3897	"O" RING, CH1000
789-094	AR-618	TRAP ASSEMBLY, MOL SIEVE, WR-12/112
789-114	AR-6317	STAINLESS PURGE TUBE FOR SC432
789-115	AR-789	LANCE TUBE ASSMBY, CS225, WO/AUTOCLNR 4.28"
789-116	AR-4174	PURGE TUBE, SC32/132
789-118	AR-9118	LANCE TUBE ASSY, 5.12LG,CS225
789-122	AR-2601	VACUUM CLEANER, SPECIFY 110/220 VAC
789-412	AR-9412	CATALYST TUBE 10.75 INCH
874-000-181	AR-4181	NICKEL BOAT, SMALL FOR LECO AMA 254 Hg

LECO® CS125/225

Description	LECO #	Alpha #
Chemical Reagents		
Cellulose	764-098	AR098
Alphacel Tungsten 2.5 kg	763-266	AR027
Alphacel Tungsten 13.6 kg	763-026	AR027B
Alphacel II Tungsten/Tin Standard Grade 2.5kg	501-008	AR008
Alphacel II Tungsten/Tin Standard Grade 13.6 kg	502-297	AR008B
Alphacel II HP - Tungsten/Tin HP 750g	5002-173	AR173
Alphacel III Tungsten 13.6 kg	763-027	AR266B
Alphacel III Tungsten 2.5kg	763-263	AR266
AlphaDri Magnesium Perchlorate 454g	501-171 501-171-HAZ	AR171
AlphaSolve A Granular 0.8 to .6mm 14 to 25 Mesh 500g	502-174 502-174-HAZ	AR2174
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Glass Wool Borosilicate Fibers 454g	501-081	AR081
Glass Wool Borosilicate Fibers 227g	763-265	AR265
Glass Wool Rope Borosilicate Fibers 454g	501-081	AR081
Glass Wool Rope Borosilicate Fibers 227g	763-265	AR265
Iron Chip Accelerator 908g	501-077	AR077
Iron Chip High Purity 454g	501-673	AR673

Iron Chip Standard Grade 11.35 kg	763-467	AR467
Platinum/Silica 15g	501-587	AR912
Tin Metal Accelerator 908g	501-076	AR076
Tin Metal Accelerator 11.35 kg	762-695	AR695
Crucibles and Boats		
Ceramic Crucible Pack of 1000	528-018	AR3818
Ceramic Crucible Pack of 500	528-050	AR3818B
Crucible Cover with 10mm Hole	528-043 619-880	AR9880
Crucible Cover with 4mm Hole Pack of 1000	528-040	AR040
Filtering Crucible Pack of 1000	528-028	AR8028
Filtering Crucible Pack of 100	528-030	AR8030
Porous Crucible Cover Pack of 1000	525-042	AR042
Reference Materials		
<u>See Alpha's Current List of Standards</u>		
O Rings and Operational Accessories		
Combustin Tube Ring Holder	772-906	AR2906
Crucible Tong	761-929	AR929
Filter Metal for manual Cleaning Head	773-311	AR311
Lance 5.2	789-118	AR9118
Lance Cleaning Rod	772-381	AR2381

Metal Scoop Assembly	773-579	AR579
O Ring	521-113	AR113
O Ring	762-058	AR058
O Ring	781-424	AR424
O Ring	772-520	AR520
O Ring Catalyst Heater Assembly	772-738	AR738
O Ring Precooler Assembly	771-991	AR991
Oscillator Tube	771-354	AR354
Pedestal Adapter	605-817	AR817
Pedestal CS125	780-953	AR922
Printer Paper Roll 4.5 x 3 x 0.375 (length 115 x diameter x internal diameter 10mm)	780-927	AR927
Printer Ribbon	780928	AR926
Secondary (screen) Filter	775-306	AR5306
Silicon Nitride Pedestal	605-816	AR816
Silicone Grease 100g	501-241	AR241
Snoop Leak Detection Fluid	502-213	AR2213
Tweezer	780-138	AR138
Wire Brush	501-062	AR022
Quartz, Glassware and Furnace Tubes		
Bubble/Filter/Reagent Tube Borosilicate Glass	771-030	AR030
Catalyst Heater Tube	772-739	AR544

Combustuion Tube	772-905	AR057
Glass Scoop	503-032	AR901

LECO® CS200/C200/S200/WC200

Description	LECO #	Alpha #
Chemical Reagents		
Cellulose	764-098	AR098
Copper Accelerator 227g	501-264	AR264
Alphacel Tungsten 2.5 kg	763-266	AR027
Alphacel Tungsten 13.6 kg	763-026	AR027B
Alphacel II Tungsten/Tin Standad Grade 2.5kg	501-008	AR008
Alphacel II Tungsten/Tin Standad Grade13.6 kg	502-297	AR008B
Alphacell II HP - Tungsten/Tin HP 750g	5002-173	AR173
Alphacel III Tungsten 13.6 kg	763-027	AR266B
Alphacel III Tungsten 2.5kg	763-263	AR266
AlphaDri Magnesium Perchlorate 454g	501-171 501-171-HAZ	AR171
AlphaSolve A Granular 0.8 to .6mm 14 to 25 Mesh 500g	502-174 502-174-HAZ	AR2174
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Glass Wool Borosilicate Fibers 454g	501-081	AR081

Glass Wool Borosilicate Fibers 227g	763-265	AR265
Glass Wool Rope Borosilicate Fibers 454g	501-081	AR081
Glass Wool Rope Borosilicate Fibers 227g	763-265	AR265
Iron Chip Accelerator 908g	501-077	AR077
Iron Chip High Purity 454g	501-673 502-230	AR673
Iron Chip Standard Grade 11.35 kg	763-467	AR467
Platinum/Silica 15g	501-587	AR912
Silica (Quartz) Wool Very Fine 50g	502-177 502-314	AR177
Tin Metal Accelerator 908g	501-076	AR076
Tin Metal Accelerator 11.35 kg	762-695	AR695
Crucibles and Boats		
Ceramic Crucible Pack of 1000	528-018	AR3818
Ceramic Crucible Pack of 500	528-050	AR3818B
Crucible Cover with 10mm Hole	528-043 619-880	AR9880
Crucible Cover with 4mm Hole Pack of 1000	528-040	AR040
Filtering Crucible Pack of 1000	528-028	AR8028
Filtering Crucible Pack of 100	528-030	AR8030
Porous Crucible Cover Pack of 1000	525-042	AR042
Reference Materials		
See Alpha's Current List of Standards		

O Rings and Operational Accessories		
Combustin Tube Ring Holder	772-906	AR2906
Crucible Tong	761-929	AR929
Filter Metal for manual Cleaning Head	773-311	AR311
Induction Coil	776-129	AR129
Lance 5.2	789-118	AR9118
Lance Cleaning Rod	772-381	AR2381
Metal Scoop Assembly	773-579	AR579
O Ring	602-920	AR920
O Ring	762-058	AR058
O Ring	781-424	AR424
O Ring	772-520	AR520
O Ring Catalyst Heater Assembly	772-738	AR738
O Ring Filter Tube	765-976	AR5967
O Ring Set	605-520-087	AR6055
Oscillator Tube	771-354	AR354
Printer Paper American Letter 8.5 x 11 x 600 sheets	601-483	AR483
Printer Ribbon	601-484	AR484
Secondary (screen) Filter	775-306	AR5306
Silicon Nitride Pedestal	605-816	AR816

Silicone Grease 100g	501-241	AR241
Snoop Leak Detection Fluid	502-213	AR2213
Tweezer	780-138	AR138
Wire Brush	501-062	AR022
Quartz, Glassware and Furnace Tubes		
Bubble/Filter/Reagent Tube Borosilicate Glass	771-030	AR030
Catalyst Heater Tube	772-739	AR544
Combustion Tube	772-905	AR057
Furnace Filter Tube Borosilicate Glass	601-390	AR390
Glass Scoop	503-032	AR901
Reagent Tube	771-373	AR373

LECO® CS230

Description	LECO #	Alpha #
Chemical Reagents		
Cellulose	764-098	AR098
Copper Accelerator 227g	501-264	AR264
Alphacel Tungsten 2.5 kg	763-266	AR027
Alphacel Tungsten 13.6 kg	763-026	AR027B
Alphacel II Tungsten/Tin Standard Grade 2.5kg	501-008	AR008
Alphacel II Tungsten/Tin Standard Grade 13.6 kg	502-297	AR008B

Alphacell II HP - Tungsten/Tin HP 750g	5002-173	AR173
Alphacel III Tungsten 13.6 kg	763-027	AR266B
Alphacel III Tungsten 2.5kg	763-263	AR266
AlphaDri Magnesium Perchlorate 454g	501-171 501-171-HAZ	AR171
AlphaSolve A Granular 0.8 to .6mm 14 to 25 Mesh 500g	502-174 502-174-HAZ	AR2174
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Glass Wool Borosilicate Fibers 454g	501-081	AR081
Glass Wool Borosilicate Fibers 227g	763-265	AR265
Glass Wool Rope Borosilicate Fibers 454g	501-081	AR081
Glass Wool Rope Borosilicate Fibers 227g	763-265	AR265
Iron Chip Accelerator 908g	501-077	AR077
Iron Chip High Purity 454g	501-673 502-230	AR673
Iron Chip Standard Grade 11.35 kg	763-467	AR467
Platinum/Silica 15g	501-587	AR912
Silica (Quartz) Wool Very Fine 50g	502-177 502-314	AR177
Crucibles and Boats		
Ceramic Crucible Pack of 1000	528-018	AR3818
Ceramic Crucible Pack of 500	528-050	AR3818B
Reference Materials		

[See Alpha's Current List of Standards](#)

O Rings and Operational Accessories

Crucible Tong	761-929	AR929
Filter Metal for manual Cleaning Head	773-311	AR311
Induction Coil	776-129	AR129
Lance 5.2	789-118	AR9118
Lance Cleaning Rod	772-381	AR2381
Metal Scoop Assembly	773-579	AR579
O Ring	602-920	AR920
O Ring	762-058	AR058
O Ring	781-424	AR424
O Ring	772-520	AR520
O Ring Catalyst Heater Assembly	772-738	AR738
O Ring Set	605-520-087	AR6055
Oscillator Tube	771-354	AR354
Pedestal Assembly	605-814	AR5814
Secondary (screen) Filter	775-306	AR5306
Silicon Nitride Pedestal	605-816	AR816
Silicone Grease 100g	501-241	AR241
Snoop Leak Detection Fluid	502-213	AR2213
Tweezer	780-138	AR138

Wire Brush	501-062	AR022
Quartz, Glassware and Furnace Tubes		
Bubble/Filter/Reagent Tube Borosilicate Glass	771-030	AR030
Catalyst Heater Tube	772-739	AR544
Combustion Tube	772-905	AR057
Furnace Filter Tube Borosilicate Glass	601-390	AR390
Glass Scoop	503-032	AR901
Reagent Tube	771-373	AR373

LECO® CS244/344/IR212/232

Description	LECO #	Alpha #
Chemical Reagents		
Cellulose	764-098	AR098
Copper Chip Accelerator 9kg	501-640	AR366
Copper Metal Accelerator 3lb	501-263	AR263
Alphacel Tungsten 2.5 kg	763-266	AR027
Alphacel Tungsten 13.6 kg	763-026	AR027B
Alphacel II Tungsten/Tin Standard Grade 2.5kg	501-008	AR008
Alphacel II Tungsten/Tin Standard Grade 13.6 kg	502-297	AR008B
Alphacel II HP - Tungsten/Tin HP 750g	5002-173	AR173
Alphacel III Tungsten 13.6 kg	763-027	AR266B

Alphacel III Tungsten 2.5kg	763-263	AR266
AlphaDri Magnesium Perchlorate 454g	501-171 501-171- HAZ	AR171
AlphaSolve A Granular 0.8 to .6mm 14 to 25 Mesh 500g	502-174 502-174- HAZ	AR2174
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Glass Wool Borosilicate Fibers 454g	501-081	AR081
Glass Wool Borosilicate Fibers 227g	763-265	AR265
Glass Wool Rope Borosilicate Fibers 454g	501-081	AR081
Glass Wool Rope Borosilicate Fibers 227g	763-265	AR265
Iron Chip Accelerator 908g	501-077	AR077
Iron Chip High Purity 454g	501-673 502-230	AR673
Iron Chip Standard Grade 11.35 kg	763-467	AR467
Platinum/Silica 15g	501-587	AR912
Silica (Quartz) Wool Very Fine 50g	502-177 502-314	AR177
Tin Metal Accelerator 908g	501-076	AR076
Tin Metal Accelerator 11.35 kg	762-695	AR695
Crucibles and Boats		
Ceramic Crucible Pack of 1000	528-018	AR3818
Ceramic Crucible Pack of 500	528-050	AR3818B

Crucible Cover with 10mm Hole	528-043 619-880	AR9880
Crucible Cover with 4mm Hole Pack of 1000	528-040	AR040
Filtering Crucible Pack of 1000	528-028	AR8028
Filtering Crucible Pack of 100	528-030	AR8030
Porous Crucible Cover Pack of 1000	525-042	AR042
Reference Materials		
See Alpha's Current List of Standards		
O Rings and Operational Accessories		
1 Gram Weight	764-078	AR4078
Brush	773-561	AR561
Combustin Tube Ring Holder	772-906	AR2906
Crucible Tong	761-929	AR929
Dust Eject Tube 19	776-095	AR095
Fan Patriot 240v	781-371	AR371
Filter Dust Box	776-234	AR6234
Filter Metal for Manual Cleaning Head	773-311	AR311
Hex Jam Nut	194-060	AR4060
Induction Coil	776-129	AR129
Lance Cleaning Rod	772-381	AR2381
Lance 5.2	789-118	AR9118

Lance Cleaning Rod	772-381	AR2381
Metal Scoop Assembly	773-579	AR579
O Ring	521-113	AR113
O Ring	578-672	AR672
O Ring	602-920	AR920
O Ring	762-058	AR058
O Ring	772-520	AR520
O Ring Catalyst Heater Assembly	772-738	AR738
O Ring Doser Valve Assembly pack of 10	772-269	AR269
O Ring Lance Tube	760-224	AR224
Oscillator Tube	771-354	AR354
Oven Fan Assembly Small	772-795	AR795
Paper Roll SC32/132/232	7700-764	AR764
Particle Filter	768-980	AR980
Pedestal Adapter	605-815	AR815
Pedestal Assembly	605-814	AR5814
Printer Ribbon CHN600/SC32/132/232/CR12	770-778	AR778
Secondary (screen) Filter	775-306	AR5306
Shield and Collar	776-922	AR6922
Silicon Nitride Pedestal	605-816	AR816
Silicone Grease 100g	501-241	AR241

Snoop Leak Detection Fluid	502-213	AR2213
Tweezer	780-138	AR138
Wire Brush	501-062	AR022
Quartz, Glassware and Furnace Tubes		
Bubble/Filter/Reagent Tube Borosilicate Glass	771-030	AR030
Catalyst Heater Tube	772-739	AR544
Ceramic Pedestal 3.05	529-007	AR529
Combustion Tube	772-905	AR057
Glass Scoop	503-032	AR901
Quartz Catalyst Tube CF10	601-064 764-720	AR064
Reagent Tube	601-065	AR1065
Reagent Tube	771-373	AR373

LECO® CS300

Description	LECO #	Alpha #
Chemical Reagents		
Cellulose	764-098	AR098
Copper Accelerator 227g	501-264	AR264
Alphacel Tungsten 2.5 kg	763-266	AR027
Alphacel Tungsten 13.6 kg	763-026	AR027B
Alphacel II Tungsten/Tin Standard Grade 2.5kg	501-008	AR008

Alphacel II Tungsten/Tin Standad Grade13.6 kg	502-297	AR008B
Alphacell II HP - Tungsten/Tin HP 750g	5002-173	AR173
Alphacel III Tungsten 13.6 kg	763-027	AR266B
Alphacel III Tungsten 2.5kg	763-263	AR266
AlphaSolve A Granular 0.8 to .6mm 14 to 25 Mesh 500g	502-174 502-174-HAZ	AR2174
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Glass Wool Borosilicate Fibers 454g	501-081	AR081
Glass Wool Borosilicate Fibers 227g	763-265	AR265
Glass Wool Rope Borosilicate Fibers 454g	501-081	AR081
Glass Wool Rope Borosilicate Fibers 227g	763-265	AR265
Iron Chip Accelerator 908g	501-077	AR077
Iron Chip High Purity 454g	501-673 502-230	AR673
Iron Chip Standard Grade 11.35 kg	763-467	AR467
Platinum/Silica 15g	501-587	AR912
Silica (Quartz) Wool Very Fine 50g	502-177 502-314	AR177
Tin Metal Accelerator 908g	501-076	AR076
Tin Metal Accelerator 11.35 kg	762-695	AR695
Crucibles and Boats		

Ceramic Crucible Pack of 1000	528-018	AR3818
Ceramic Crucible Pack of 500	528-050	AR3818B
Crucible Cover with 10mm Hole	528-043 619-880	AR9880
Crucible Cover with 4mm Hole Pack of 1000	528-040	AR040
Filtering Crucible Pack of 1000	528-028	AR8028
Filtering Crucible Pack of 100	528-030	AR8030
Porous Crucible Cover Pack of 1000	525-042	AR042
Reference Materials		
See Alpha's Current List of Standards		
O Rings and Operational Accessories		
Combustion Tube Ring Holder	772-906	AR2906
Crucible Tong	761-929	AR929
Filter Metal for Manual Cleaning Head	773-311	AR311
Induction Coil	776-129	AR129
Lance 5.2	789-118	AR9118
Lance Cleaning Rod	772-381	AR2381
Metal Scoop Assembly	773-579	AR579
O Ring	602-920	AR920
O Ring	762-058	AR058

O Ring	781-424	AR424
O Ring	772-520	AR520
O Ring Catalyst Heater Assembly	772-738	AR738
O Ring Filter Tube	765-976	AR5967
O Ring Set	605-520-087	AR6055
Oscillator Tube	771-354	AR354
Printer Paper American Letter 8.5 x 11 x 600 sheets	601-483	AR483
Printer Ribbon	601-484	AR484
Secondary (screen) Filter	775-306	AR5306
Silicon Nitride Pedestal	605-816	AR816
Silicone Grease 100g	501-241	AR241
Snoop Leak Detection Fluid	502-213	AR2213
Tweezer	780-138	AR138
Wire Brush	501-062	AR022
Quartz, Glassware and Furnace Tubes		
Bubble/Filter/Reagent Tube Borosilicate Glass	771-030	AR030
Catalyst Heater Tube	772-739	AR544
Combustion Tube	772-905	AR057
Furnace Filter Tube Borosilicate Glass	601-390	AR390
Glass Scoop	503-032	AR901
Reagent Tube	771-373	AR373

LECO® CS400/C400/S400

Description	LECO #	Alpha #
Chemical Reagents		
Cellulose	764-098	AR098
Copper Accelerator 227g	501-264	AR264
Alphacel Tungsten 2.5 kg	763-266	AR027
Alphacel Tungsten 13.6 kg	763-026	AR027B
Alphacel II Tungsten/Tin Standard Grade 2.5kg	501-008	AR008
Alphacel II Tungsten/Tin Standard Grade 13.6 kg	502-297	AR008B
Alphacel II HP - Tungsten/Tin HP 750g	5002-173	AR173
Alphacel III Tungsten 13.6 kg	763-027	AR266B
Alphacel III Tungsten 2.5kg	763-263	AR266
AlphaDri Magnesium Perchlorate 454g	501-171 501-171-HAZ	AR171
AlphaSolve A Granular 0.8 to .6mm 14 to 25 Mesh 500g	502-174 502-174-HAZ	AR2174
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Glass Wool Borosilicate Fibers 454g	501-081	AR081
Glass Wool Borosilicate Fibers 227g	763-265	AR265
Glass Wool Rope Borosilicate Fibers 454g	501-081	AR081
Glass Wool Rope Borosilicate Fibers 227g	763-265	AR265

Iron Chip Accelerator 908g	501-077	AR077
Iron Chip High Purity 454g	501-673 502-230	AR673
Iron Chip Standard Grade 11.35 kg	763-467	AR467
Platinum/Silica 15g	501-587	AR912
Silica (Quartz) Wool Very Fine 50g	502-177 502-314	AR177
Tin Metal Accelerator 908g	501-076	AR076
Tin Metal Accelerator 11.35 kg	762-695	AR695
Crucibles and Boats		
Ceramic Crucible Pack of 1000	528-018	AR3818
Ceramic Crucible Pack of 500	528-050	AR3818B
Crucible Cover with 10mm Hole	528-043 619-880	AR9880
Crucible Cover with 4mm Hole Pack of 1000	528-040	AR040
Filtering Crucible Pack of 1000	528-028	AR8028
Filtering Crucible Pack of 100	528-030	AR8030
Porous Crucible Cover Pack of 1000	525-042	AR042
Reference Materials		
See Alpha's Current List of Standards		
O Rings and Operational Accessories		
1 Gram Weight	764-078	AR4078
Combustion Tube Ring Holder	772-906	AR2906

Crucible Tong	761-929	AR929
Filter Dust Box	776-234	AR6234
Filter Metal for Manual Cleaning Head	773-311	AR311
Hex Jam Nut	194-060	AR4060
Induction Coil	776-129	AR129
Lance 5.2	789-118	AR9118
Lance Cleaning Rod	772-381	AR2381
Metal Scoop Assembly	773-579	AR579
O Ring	602-920	AR920
O Ring	762-058	AR058
O Ring	781-424	AR424
O Ring	772-520	AR520
O Ring Catalyst Heater Assembly	772-738	AR738
O Ring Filter Tube	765-976	AR5967
O Ring Precooler Assembly	771-991	AR991
O Ring Set	605-520-087	AR6055
Oscillator Tube	771-354	AR354
Printer Paper American Letter 8.5 x 11 x 600 sheets	601-483	AR483
Printer Ribbon	601-484	AR484
Secondary (screen) Filter	775-306	AR5306

Silicon Nitride Pedestal	605-816	AR816
Silicone Grease 100g	501-241	AR241
Snoop Leak Detection Fluid	502-213	AR2213
Tweezer	780-138	AR138
Wire Brush	501-062	AR022
Quartz, Glassware and Furnace Tubes		
Bubble/Filter/Reagent Tube Borosilicate Glass	771-030	AR030
Catalyst Heater Tube	772-739	AR544
Combustion Tube	772-905	AR057
Furnace Filter Tube Borosilicate Glass	601-390	AR390
Glass Scoop	503-032	AR901
Reagent Tube	771-373	AR373

LECO® CS444/IR4112/432

Description	LECO #	Alpha #
Chemical Reagents		
Cellulose	764-098	AR098
Copper Chip Accelerator 9kg	501-640	AR366
Copper Metal Accelerator 3lb	501-263	AR263
Alphacel Tungsten 2.5 kg	763-266	AR027
Alphacel Tungsten 13.6 kg	763-026	AR027B

Alphacel II Tungsten/Tin Standad Grade 2.5kg	501-008	AR008
Alphacel II Tungsten/Tin Standad Grade13.6 kg	502-297	AR008B
Alphacell II HP - Tungsten/Tin HP 750g	5002-173	AR173
Alphacel III Tungsten 13.6 kg	763-027	AR266B
Alphacel III Tungsten 2.5kg	763-263	AR266
AlphaDri Magnesium Perchlorate 454g	501-171 501-171-HAZ	AR171
AlphaSolve A Granular 0.8 to .6mm 14 to 25 Mesh 500g	502-174 502-174-HAZ	AR2174
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Glass Wool Borosilicate Fibers 454g	501-081	AR081
Glass Wool Borosilicate Fibers 227g	763-265	AR265
Glass Wool Rope Borosilicate Fibers 454g	501-081	AR081
Glass Wool Rope Borosilicate Fibers 227g	763-265	AR265
Iron Chip Accelerator 908g	501-077	AR077
Iron Chip High Purity 454g	501-673 502-230	AR673
Iron Chip Standard Grade 11.35 kg	763-467	AR467
Platinum/Silica 15g	501-587	AR912
Rare Earth Copper Oxide Granular 0.85 to 1.7mm 50g	501-170	AR170
Silica (Quartz) Wool Very Fine 50g	502-177 502-314	AR177
Tin Metal Accelerator 908g	501-076	AR076

Tin Metal Accelerator 11.35 kg	762-695	AR695
Crucibles and Boats		
Ceramic Crucible Pack of 1000	528-018	AR3818
Crucible Cover with 10mm Hole	528-043 619-880	AR9880
Crucible Cover with 4mm Hole Pack of 1000	528-040	AR040
Filtering Crucible Pack of 1000	528-028	AR8028
Filtering Crucible Pack of 100	528-030	AR8030
Porous Crucible Cover Pack of 1000	525-042	AR042
Reference Materials		
See Alpha's Current List of Standards		
O Rings and Operational Accessories		
1 Gram Weight	764-078	AR4078
Brush	604-163	AR163
Brush	773-314	AR3314
Brush	604-271	AR271
Combustion Tube Ring Holder	772-906	AR2906
Crucible Tong	761-929	AR929
Filter Metal for Autoclean Head	778-596	AR596
Hex Jam Nut	194-060	AR4060
Induction Coil	776-129	AR129

Lance 4.25	789-115	AR789
Lance Cleaning Rod	772-381	AR2381
Metal Scoop Assembly	773-579	AR579
O Ring	602-920	AR920
O Ring	604-172	AR172
O Ring	762-058	AR058
O Ring	772-520	AR520
O Ring Catalyst Heater Assembly	772-738	AR738
O Ring Doser Valve Assembly pack of 10	772-269	AR269
O Ring Precooler Assembly	771-991	AR991
O Ring Reagent Tube	773-141 773-320	AR320
Particle Filter	768-980	AR980
Pedestal Adapter	605-815	AR815
Pedestal Assembly	605-814	AR5814
Printer Paper American Letter 8.5 x 11 x 600 sheets	601-483	AR483
Printer Ribbon	601-484	AR484
Secondary (screen) Filter	775-306	AR5306
Silicon Nitride Pedestal	605-816	AR816
Silicone Grease 100g	501-241	AR241
Snoop Leak Detection Fluid	502-213	AR2213
Tweezer	780-138	AR138

Vacuum Filter CS444 IR431 & IR412	782-880	AR2880
Wire Brush	501-062	AR022
Quartz, Glassware and Furnace Tubes		
Bubble/Filter/Reagent Tube Borosilicate Glass	771-030	AR030
Catalyst Heater Tube	772-739	AR544
Combustion Tube	772-905	AR057
Glass Scoop	503-032	AR901
Quartz Catalyst Tube	604-112	AR604
Reagent Tube	761-047	AR047
Reagent Tube	771-373	AR373

LECO® CS600

Description	LECO #	Alpha #
Chemical Reagents		
Cellulose	764-098	AR098
Copper Chip Accelerator 9kg	501-640	AR366
Copper Metal Accelerator 3lb	501-263	AR263
Alphacel Tungsten 2.5 kg	763-266	AR027
Alphacel Tungsten 13.6 kg	763-026	AR027B
Alphacel II Tungsten/Tin Standard Grade 2.5kg	501-008	AR008

Alphacel II Tungsten/Tin Standad Grade13.6 kg	502-297	AR008B
Alphacell II HP - Tungsten/Tin HP 750g	5002-173	AR173
Alphacel III Tungsten 13.6 kg	763-027	AR266B
Alphacel III Tungsten 2.5kg	763-263	AR266
AlphaDri Magnesium Perchlorate 454g	501-171 501-171- HAZ	AR171
AlphaSolve A Granular 0.8 to .6mm 14 to 25 Mesh 500g	502-174 502-174- HAZ	AR2174
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Glass Wool Borosilicate Fibers 454g	501-081	AR081
Glass Wool Borosilicate Fibers 227g	763-265	AR265
Glass Wool Rope Borosilicate Fibers 454g	501-081	AR081
Glass Wool Rope Borosilicate Fibers 227g	763-265	AR265
Iron Chip Accelerator 908g	501-077	AR077
Iron Chip High Purity 454g	501-673 502-230	AR673
Iron Chip Standard Grade 11.35 kg	763-467	AR467
Platinum/Silica 15g	501-587	AR912
Rare Earth Copper Oxide Granular 0.85 to 1.7mm 50g	501-170	AR170
Silica (Quartz) Wool Very Fine 50g	502-177 502-314	AR177
Tin Metal Accelerator 908g	501-076	AR076
Tin Metal Accelerator 11.35 kg	762-695	AR695

Crucibles and Boats		
Ceramic Crucible Pack of 1000	528-018	AR3818
Ceramic Crucible Pack of 500	528-050	AR3818B
Crucible Cover with 10mm Hole	528-043 619-880	AR9880
Crucible Cover with 4mm Hole Pack of 1000	528-040	AR040
Filtering Crucible Pack of 1000	528-028	AR8028
Filtering Crucible Pack of 100	528-030	AR8030
Porous Crucible Cover Pack of 1000	525-042	AR042
Reference Materials		
See Alpha's Current List of Standards		
O Rings and Operational Accessories		
Brush	604-163	AR163
Brush	604-271	AR271
Combustion Tube Ring Holder	772-906	AR2906
Crucible Tong	761-929	AR929
Filter Metal for Autoclean Head	778-596	AR596
Hex Jam Nut	194-060	AR4060
Induction Coil	776-129	AR129
Lance 4.25	789-115	AR789

Lance Cleaning Rod	772-381	AR2381
Metal Scoop Assembly	773-579	AR579
O Ring	602-920	AR920
O Ring	604-172	AR172
O Ring	617-083	AR7083
O Ring	762-058	AR058
O Ring	772-520	AR520
O Ring Catalyst Heater Assembly	772-738	AR738
O Ring Doser Valve Assembly pack of 10	772-269	AR269
O Ring Filter Tube	765-976	AR5967
O Ring Precooler Assembly	771-991	AR991
O Ring Reagent Tube	773-141 773-320	AR320
Oscillator Tube	771-354	AR354
Particle Filter	768-980	AR980
Pedestal Adapter	605-815	AR815
Pedestal Assembly	605-814	AR5814
Printer Paper American Letter 8.5 x 11 x 600 sheets	601-483	AR483
Printer Ribbon	601-484	AR484
Screen Mesh SC32/132/232	760-193	AR193
Secondary (screen) Filter	775-306	AR5306
Silicon Nitride Pedestal	605-816	AR816

Silicone Grease 100g	501-241	AR241
Snoop Leak Detection Fluid	502-213	AR2213
Tweezer	780-138	AR138
Vacuum Filter CS444 IR431 & IR412	782-880	AR2880
Wire Brush	501-062	AR022
Quartz, Glassware and Furnace Tubes		
Bubble/Filter/Reagent Tube Borosilicate Glass	771-030	AR030
Catalyst Heater Tube	772-739	AR544
Combustion Tube	772-905	AR057
Filling Funnel Borosilicate Glass	502-023	AR473
Furnace Filter Tube Borosilicate Glass	601-390	AR390
Glass Scoop	503-032	AR901
Quartz Catalyst Tube	604-112	AR604
Reagent Tube	771-373	AR373

LECO® DH603

Description	LECO #	Alpha #
Chemical Reagents		
Copper Oxide Fine Wires 4 x 0.5mm 100g	502-190	AEB1002

AlphaDri Magnesium Perchlorate 454g	501-171 501-171- HAZ	AR171
AlphaSolve A Granular 0.8 to .6mm 14 to 25 Mesh 500g	502-174 502-174-HAZ	AR2174
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Glass Wool Borosilicate Fibers 454g	501-081	AR081
Glass Wool Borosilicate Fibers 227g	763-265	AR265
Glass Wool Rope Borosilicate Fibers 454g	501-081	AR081
Glass Wool Rope Borosilicate Fibers 227g	763-265	AR265
Quartz Wool Wadding 200g	501-608	AR1608
Reference Materials		
See Alpha's Current List of Standards		
O Rings and Operational Accessories		
O Ring	617-083	AR7083
O Ring	772-520	AR520
O Ring Filter Tube	765-976	AR5967
Secondary (screen) Filter	775-306	AR5306
Silicone Grease 100g	501-241	AR241
Snoop Leak Detection Fluid	502-213	AR2213
Teflon Disconnect block	773-931	AR9931
Tweezer	780-138	AR138
Quartz, Glassware and Furnace Tubes		

Catalyst Tube	789-412	AR9412
Reagent Tube	617-548	AR7548

LECO® RH402

Description	LECO #	Alpha #
Tin Capsules for Liquids		
Smooth Wall Tin Capsules Round Base 18 x 5mm pack of 100	501-089	AR059
Smooth Wall Tin Capsules Round Base 18 x 6mm pack of 100	502-040	AR2166
Chemical Reagents		
AlphaDri Magnesium Perchlorate 454g	501-171 501-171- HAZ	AR171
AlphaSolve A Granular 0.8 to .6mm 14 to 25 Mesh 500g	502-174 502-174- HAZ	AR2174
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Glass Wool Borosilicate Fibers 454g	501-081	AR081
Glass Wool Borosilicate Fibers 227g	763-265	AR265
Glass Wool Rope Borosilicate Fibers 454g	501-081	AR081
Glass Wool Rope Borosilicate Fibers 227g	763-265	AR265
Schutze Reagent 100g	761-747	AR611

Crucibles and Boats		
Quartz Crucible pack of 24	603-710	AR710
Short Graphite Crucible Degassed pack of 100	769-520	AR9520C
Reference Materials		
See Alpha's Current List of Standards		
O Rings and Operational Accessories		
Air Filter Coating	766-646	AR6646
Crucible Tweezers	766-053	AR766
Fan Patriot 240v	781-371	AR371
Knurled Ring Collar	774-690	AR690
O Ring	602-803	AR803
O Ring	762-357	AR357
O Ring	602-802	AR802
O Ring	603-910	AR3910
O Ring	772-520	AR520
O Ring Catalyst Heater Assembly	772-738	AR738
O Ring Catalyst Heater Tube	778-114	AR8114
O Ring Sample Drop	772-910	AR2910
Particle Filter	768-980	AR980
Preparation File	501-234	AR234
Printer Paper American Letter 8.5 x 11 x 600 sheets	601-483	AR483

Printer Ribbon	601-484	AR484
Quad O Ring	603-970	AR3970
Secondary (screen) Filter	775-306	AR5306
Silicone Grease 100g	501-241	AR241
Snoop Leak Detection Fluid	502-213	AR2213
Tweezer	780-138	AR138
Up/Down Solenoid	776-020 776-685	AR6685
Quartz, Glassware and Furnace Tubes		
Catalyst Heater Tube Transparent Silica	776-279	AR279
Funnel	767-473	AR473
Reaction Tube	604-241	AR4241
Reaction Tube RH402	603-707	AR707
Reagent Tube Borosilicate Glass	775-601	AR601

LECO® RH404

Description	LECO #	Alpha #
Tin Capsules for Liquids		
Smooth Wall Tin Capsules Round Base 18 x 5mm pack of 100	501-089	AR059
Smooth Wall Tin Capsules Round Base 18 x 6mm pack of 100	502-040	AR2166

Chemical Reagents		
Copper Turnings Fine Long Strands Reduced and Annealed 200g	501-621 502-258 502-295	AR621
AlphaDri Magnesium Perchlorate 454g	501-171 501-171- HAZ	AR171
AlphaSolve A Granular 0.8 to .6mm 14 to 25 Mesh 500g	502-174 502-174-HAZ	AR2174
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Glass Wool Borosilicate Fibers 454g	501-081	AR081
Glass Wool Borosilicate Fibers 227g	763-265	AR265
Glass Wool Rope Borosilicate Fibers 454g	501-081	AR081
Glass Wool Rope Borosilicate Fibers 227g	763-265	AR265
Crucibles and Boats		
Graphite Crucible pack of 100	764-330	AR331S
High Temp Graphite Crucible pack of 100	782-720	AR720S
High Temp Graphite Crucible pack of 1000	782-720	AR782HD
Reference Materials		
See Alpha's Current List of Standards		
O Rings and Operational Accessories		
1 Gram Weight	764-078	AR4078
Air Filter Coating	766-646	AR6646

Brush	782-938	AR938
Brush Handle		AR938H
Filter Metal for Manual Cleaning Head	773-311	AR311
Brush Head	77-652	AR652
Cleaning Brush Assembly	774-649	AR653
Cleaning Brush Assembly	782-918	AR877H
Cleaning Brush Assembly	782-942	AR938H
Crucible Tweezers	760-139	AR766
Crucible Tweezers	766-053	AR766
Knurled Ring Collar	774-690	AR690
Lower Electrode Tip	782-721	AR721
Lubricool Cooling Agent	765-933	AR5933
O Ring	568-016 760-226	AR428
O Ring	771-407	AR407
O Ring	775-312	AR312
O Ring	772-520	AR520
O Ring Catalyst Heater Assembly	772-738	AR738
O Ring Catalyst Heater Tube	778-114	AR8114
O Ring Combustion System	774-489	AR489
Particle Filter	768-980	AR980

Preparation File	501-234	AR234
Printer Paper American Letter 8.5 x 11 x 600 sheets	601-483	AR483
Printer Ribbon	601-484	AR484
Secondary (screen) Filter	775-306	AR5306
Silicone Grease 100g	501-241	AR241
Small Cleaning Brush Assembly	775-416	AR293H
Snoop Leak Detection Fluid	502-213	AR2213
Trim Step Brush	782-877	AR877
Tweezer	780-138	AR138
Up/Down Solenoid	776-020 776-685	AR6685
Vacuum Filter CS444 IR431 & IR412	782-880	AR2880
Quartz, Glassware and Furnace Tubes		
Catalyst Heater Tube Transparent Silica	776-279	AR279
Funnel	767-473	AR473
Glass Scoop	503-032	AR901
Reagent Tube	761-047	AR047
Reagent Tube Borosilicate Glass	775-601	AR601

LECO® RO478

Description	LECO #	Alpha #
Tin Capsules/Pans for Solids		
Tin Capsules Pressed Standard Weight 4 x 3.2mm pack of 500	502-227	ATD1000
Chemical Reagents		
Copper Turnings Fine Long Strands Reduced and Annealed 200g	501-621 502-258 502-295	AR621
AlphaDri Magnesium Perchlorate 454g	501-171 501-171-HAZ	AR171
AlphaSolve A Granular 0.8 to .6mm 14 to 25 Mesh 500g	502-174 502-174-HAZ	AR2174
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Glass Wool Borosilicate Fibers 454g	501-081	AR081
Glass Wool Borosilicate Fibers 227g	763-265	AR265
Glass Wool Rope Borosilicate Fibers 454g	501-081	AR081
Glass Wool Rope Borosilicate Fibers 227g	763-265	AR265
Crucibles and Boats		
High Temp Graphite Crucible pack of 1000	782-720-S	AR720S
High Temp Graphite Crucible pack of 1000	782-720-S	AR720S
High Temp Graphite Crucible pack of 5000	782-720-S	AR720S
Ultra High Purity Nickel Basket 1g pack of 100	502-344	AR2344
O Rings and Operational Accessories		

O Ring	772-520	AR520
O Ring Catalyst Heater Assembly	772-738	AR738
O Ring Catalyst Heater Tube	778-114	AR8114
O Ring Sample Drop Block	589-551	AR422
O Ring Sample Drop	772-910	AR2910
O Ring Sample Drop Block Piston	774-613	AR4613
Printer Paper American Letter 8.5 x 11 x 600 sheets	601-483	AR483
Printer Ribbon	601-484	AR484
Secondary (screen) Filter	775-306	AR5306
Silicone Grease 100g	501-241	AR241
Snoop Leak Detection Fluid	502-213	AR2213
Quartz, Glassware and Furnace Tubes		
Catalyst Heater Tube Transparent Silica	776-279	AR279
Filling Funnel Borosilicate Glass	502-023	AR473
Reagent Tube Borosilicate Glass	775-601	AR601

LECO® TC436/TN414/RO416

Description	LECO #	Alpha #
Nickel Capsules for Liquids		
Nickel Capsules 11 x 5mm pack of 100	502-183	AR2183
Nickel Capsules 6.5 x 10mm pack of 100	502-183	AR2183
Tin Capsules for Liquids		
Smooth Wall Tin Capsules Round Base 18 x 5mm pack of 100	501-089	AR059
Smooth Wall Tin Capsules Round Base 18 x 6mm pack of 100	502-040	AR2166
Chemical Reagents		
Copper Turnings Fine Long Strands Reduced and Annealed 200g	501-621 502-258 502-295	AR621
Copper Wires Fine Wires Reduced 4 x 0.5mm 100g	502-189	AR189
Copper wires Pound Pack Fine Wires Reduced 4 x 0.5mm 454g	502-304-500	AR2304-500
AlphaDri Magnesium Perchlorate 454g	501-171 501-171-HAZ	AR171
AlphaSolve A Granular 0.8 to .6mm 14 to 25 Mesh 500g	502-174 502-174-HAZ	AR2174
Alphasolve B Granular 1.6 to 3mm 7 to 14 Mesh 500g	502-176	AR2176
Glass Wool Borosilicate Fibers 454g	501-081	AR081
Glass Wool Borosilicate Fibers 227g	763-265	AR265
Glass Wool Rope Borosilicate Fibers 454g	501-081	AR081

Glass Wool Rope Borosilicate Fibers 227g	763-265	AR265
Graphite Powder 100g	501-073	AR073
Nibbled Nickel Flux 100g	501-598	AR598
Rare Earth Copper Oxide Granular 0.85 to 1.7mm 50g	501-170	AR170
Crucibles and Boats		
Graphite Crucible pack of 100	776-247	AR053P
Graphite Crucible pack of 1000	776-247	AR6247
High Purity Inner Crucible pack of 1000	785-795HD	AR2795
High Purity Outer Crucible pack of 100	782-703-100	AR2703
High Temp Graphite Crucible pack of 100	782-720	AR720S
High Temp Graphite Crucible pack of 1000	782-720	AR782HD
High Temp Graphite Crucible pack of 100	782-720-S	AR720S
High Temp Graphite Crucible pack of 1000	782-720-S	AR720S
High Temp Graphite Crucible pack of 5000	782-720-S	AR720S
Inner Graphite Crucible pack of 100	775-431	AR5431
Inner Graphite Crucible pack of 100	775-892	AR431P
Nickel Basket 1.5g Each pack of 100	763-029	AR2345
Nickel Basket 1g Each pack of 100	763-065	AR2344
Outer Graphite Crucible pack of 50	775-433	AR433
Ultra High Purity Nickel Basket 1g pack of 100	502-344	AR2344

Reference Materials

[See Alpha's Current List of Standards](#)

O Rings and Operational Accessories

1 Gram Weight	764-078	AR4078
Brush	782-938	AR938
Brush Handle		AR938H
Brush Head	188-012	AR293
Brush Head	77-652	AR652
Cleaning Brush Assembly	774-649	AR653
Cleaning Brush Assembly	782-918	AR877H
Cleaning Brush Assembly	782-942	AR938H
Crucible Tweezers	766-053	AR766
Knurled Ring Collar	774-690	AR690
Lower Electrode Tip	618-482 620-814 774-691	AR691
Lower Electrode Tip	782-721	AR721
Lubricool Cooling Agent	765-933	AR5933
O Ring	568-016 760-226	AR428
O Ring	771-407	AR407
O Ring	775-312	AR312
O Ring	772-520	AR520

O Ring Catalyst Heater Assembly	772-738	AR738
O Ring Catalyst Heater Tube	778-114	AR8114
O Ring Combustion System	774-489	AR489
O Ring Sample Drop Block Piston	774-613	AR4613
Particle Filter	768-980	AR980
Preparation File	501-234	AR234
Printer Paper American Letter 8.5 x 11 x 600 sheets	601-483	AR483
Printer Ribbon	601-484	AR484
Secondary (screen) Filter	775-306	AR5306
Silicone Grease 100g	501-241	AR241
Small Cleaning Brush Assembly	775-416	AR293H
Snoop Leak Detection Fluid	502-213	AR2213
Trim Step Brush	782-877	AR877
Tweezer	780-138	AR138
Vacuum Filter CS444 IR431 & IR412	782-880	AR2880
Quartz, Glassware and Furnace Tubes		
Catalyst Heater Tube	775-355	AR355
Catalyst Heater Tube Transparent Silica	776-279	AR279
Filling Funnel Borosilicate Glass	502-023	AR473
Funnel	767-473	AR473

Glass Scoop	503-032	AR901
Reagent Tube	761-047	AR047
Reagent Tube	771-373	AR373
Reagent Tube Borosilicate Glass	775-601	AR601