



Consommables & appareils de préparation des échantillons pour analyses par Fluorescence de rayons X

Choisissez votre Gamme

[Coupelles PP pour liquide et poudre](#)

[Films pour coupelles](#)

[Coupelles Aluminium pour pastillage](#)

[Liants pour pastillage](#)

[Presse à pastiller](#)

[Moules de compression](#)

[Machine de fusion gaz et électrique
pour la fusion de perles](#)

[Fondants pour la fusion des perles](#)

[Verres de calibration](#)

Vous serez automatiquement redirigé à la section souhaitée

Couppelles à échantillons montées sans film

SC-3332	Sample cup 32mm	Sample cup cap designed for Spectro XEPOS X-LAB 2000
SC-3332-240	Sample cup 32mm + film Polypropylene 4,0µm	Sample cup designed for Spectro XEPOS and X-LAB with cap outerdiam 32mm already mounted with Polypropylene 4µm (TF-240) set=100pcs
SC-3335	Sample cup 35mm	Cup ring cap for Philips/Panalytical
SC-3340	Sample cup 40mm	Inner ring outer ring cap for Bruker SRS S4 S8 ARL Optimix
SC-3340T	Sample cup 40mm	Inner ring outer ring cap for Bruker SRS S4 S8
SC-3345	Sample cup 45mm	Inner ring outer ring cap universal for 45mm
SC-3351	Sample cup 51mm	Cup ring cap universal for 51.5mm holder
SC-4031	Sample cup 31mm	Cup ring single open ended universal open for 32mm holder
SC-4047	Sample cup 47mm	Cup ring single open ended for Horiba SOLD OUT!
SC-4131	Sample cup 31mm	Cup ring collar single open ended universal for 32mm holder
SC-4140	Sample cup 40mm	Cup ring collar single open ended universal for 40mm holder
SC-4231	Sample cup 32mm	Cup 2 rings double open ended universal for 32mm holder
SC-4232	Sample cup 32mm	Cup 2 rings double open ended for Spectro XEPOS X-LAB 2000
SC-4240	Sample cup 40mm	Cup 2 rings double open ended universal for 40mm holder
SC-4331	Sample cup 31mm	Cup ring cap universal for 32mm holder
SC-4332	Sample cup 32mm	Cup ring cap for Spectro XEPOS X-LAB 2000
SC-4335	Sample cup 35mm	Cup ring cap for Philips/Panalytical
SC-4340	Sample cup 40mm	Cup ring cap universal for 40mm holder

SC-4345	Sample cup 45mm	Cup ring cap universal for 45mm holder
SC-7332	Sample cup 34mm	Inner cup and cap designed for Oxford
SC-7332cap	Sample cup cap 32mm	Only cap designed for Oxford LABX
SC-7332inner	Sample cup Inner 32mm	Only inner cup designed for Oxford LABX
SC-8047	Sample cup 47mm	Cup ring single open ended for Horiba 47mm holder
SC-MC-1520	Sample cup 31mm	XRF micro sample cups double open ended with vented caps

Coupelles à échantillons montées avec film

SC-3332-125	Sample cup 32mm + film Mylar 2,5 μ m	Sample cup designed for Spectro XEPOS and X-LAB with cap outerdiam 32mm already mounted with Mylar 2,5 μ m (TF-125) set=100pcs SOLD OUT!
SC-3332-160	Sample cup 32mm + film Mylar 6,0 μ m	Sample cup designed for Spectro XEPOS and X-LAB with cap outerdiam 32mm already mounted with Mylar 6,0 μ m (TF-160) set=100pcs
SC-3332-260	Sample cup 32mm + film Polypropylene 6,0 μ m	Sample cup designed for Spectro XEPOS and X-LAB with cap outerdiam 32mm already mounted with polypropylene 6,0 μ m (TF-260) set=100pcs
SC-3335-160	Sample cup 35mm + film Mylar 6,0 μ m	Sample cup designed for Philips/Panalytical with cap outerdiam 35mm already mounted with mylar 6,0 μ m (TF-160) set=100pcs
SC-3335-240	Sample cup 35mm + film Polypropylene 4,0 μ m	Sample cup designed for Philips/Panalytical with cap outerdiam 35mm already mounted with polypropylene 4,0 μ m (TF-240) set=100pcs
SC-3335-260	Sample cup 35mm + film Polypropylene 6,0 μ m	Sample cup designed for Philips/Panalytical with cap outerdiam 35mm already mounted with polypropylene 6,0 μ m (TF-260) set=100pcs
SC-3340-240	Sample cup 40mm + film Polypropylene 4,0 μ m	Sample cup designed for Bruker SRS,S4,S8,ARL Optimix with cap outerdiam 35mm already mounted with polypropylene 4,0 μ m (TF-240) set=100pcs
SC-	Sample cup 40mm	Sample cup designed for Bruker, SRS,S4,S8 ARL Optimix with cap

3340-260	+ film Polypropylene 6,0µm	outerdiam 40mm already mounted with polypropylene 6,0µm (TF-260) set=100pcs
SC-3340T-240	Sample cup 40mm + film Polypropylene 4,0µm	Sample cup designed for Bruker SRS,S4,S8,ARL Optimix with cap outerdiam 40mm already mounted with polypropylene 4,0µm (TF-240) set=100pcs
SC-3345-160	Sample cup 45mm + film Polypropylene 6,0µm	Sample cup with cap universal outerdiam 45mm already mounted with polypropylene 6,0µm (TF-160) set=100pcs
SC-3345-240	Sample cup 45mm + film Polypropylene 4,0µm	Sample cup with cap universal outerdiam 45mm already mounted with polypropylene 4,0µm (TF-240) set=100pcs
SC-4331-135	Sample cup 32mm + film Mylar 3,5µm	Sample cup with cap universal outerdiam 32mm already mounted with mylar 3,5µm (TF-135) set=100pcs
SC-4331-160	Sample cup 32mm + film Mylar 6,0µm	Sample cup with cap universal outerdiam 32mm already mounted with mylar 6,0µm (TF-160) set=100pcs
SC-4331-240	Sample cup 32mm + film Polypropylene 4,0µm	Sample cup with cap universal outerdiam 32mm already mounted with polypropylene 6,0µm (TF-240) set=100pcs
SC-4331-260	Sample cup 32mm + film Polypropylene 6,0µm	Sample cup with cap universal outerdiam 32mm already mounted with polypropylene 6,0µm (TF-260) set=100pcs
SC-4332-160	Sample cup 32mm + film Mylar 6,0µm	Sample cup designed for Spectro XEPOS and X-LAB with cap outerdiam 32mm already mounted with mylar 6,0µm (TF-160) set=100pcs
SC-4332-240-115	Sample cup 32mm + film Mylar	Sample cup designed for Spectro XEPOS and X-LAB with cap outerdiam 32mm already mounted with 2 films Mylar 4,0µm (TF-240) and Mylar 1,5µm (TF115) set=100pcs
SC-4332-260	Sample cup 32mm + film Polypropylene 6,0µm	Sample cup designed for Spectro XEPOS and X-LAB with cap outerdiam 32mm already mounted with polypropylene 6,0µm (TF-260) set=100pcs
SC-4340-240	Sample cup 40mm + film Polypropylene 4,0µm	Sample cup with cap universal outerdiam 40mm already mounted with polypropylene 4,0µm (TF-240) set=100pcs

Films XRF

12345	Thin Film Microporous Polypropylene 6.4cm x 5.1m Roll	Microporous Polypropylene 6.4cm x 5.1m Roll
TF-025	Thin Film Microporous Polypropylene 6.4cm x 5.1m Roll	Microporous Polypropylene 6.4cm x 5.1m Roll
TF-050	Thin Film Microporous Teflon 6.4cm x 5m Roll	Microporous Teflon 6.4cm x 5,1m Roll
TF-112	Thin Film Mylar 12 μ m 7.6cm x 91.4m Roll	Mylar 12 μ m 7.6cm x 91.4m Roll
TF-115	Thin Film Mylar 1.5 μ m 7.6cm x 91.4m Roll	Mylar 1.5 μ m 7.6cm x 91.4m Roll
TF-115-345	Thin Film Mylar 1.5 μ m 7.6cm x 7.6cm Square	Mylar 1.5 μ m Cut Sheets 500 Pieces
TF-125	Thin Film Mylar 2.5 μ m 7.6cm x 91.4m Roll	Mylar 2.5 μ m 7.6cm x 91.4m Roll
TF-125-255	Thin Film Mylar 2.5 μ m 6.4cm Circle	Mylar 2.5 μ m 6.4cm 1Set=500 Pieces
TF-125-345	Thin Film Mylar 2.5 μ m 7.6cmx7.6cm Square	Mylar 2.5 μ m Cut Sheets 500 Pieces
TF-1335	Thin Film Hostaphan -NO LONGER AVAILABLE-	ALTERNATIVE: TF-135
TF-1335-345	Thin Film Hostaphan -NO LONGER AVAILABLE-	ALTERNATIVE: TF-135-345
TF-135	Thin Film Mylar 3.5 μ m 7.6cm x 91.4m Roll	Mylar 3.5 μ m 7.6cm x 91.4m Roll
TF-135-255	Thin Film Mylar 3.5 μ m D6.4 Circle	Mylar 3.5 μ m D6.4cm 1Set=500Pieces
TF-135-345	Thin Film Mylar 3.5 μ m 7.6cm x 7.6cm Square	Mylar 3.5 μ m Cut Sheets 1Set=500 Pieces
TF-135-355	Thin Film Mylar 3.5 μ m D8.9cm Circle	Mylar 3.5 μ m D8.9cm 1Set=500 Pieces
TF-160	Thin Film Mylar 6.0 μ m 7.6cm x 91.4m Roll	Mylar 6.0 μ m 7.6cm x 91.4m Roll
TF-160-255	Thin Film Mylar 6.0 μ m D6.4cm Circle	Mylar 6.0 μ m D6.4cm 1Set=500 Pieces
TF-160-345	Thin Film Mylar 6.0 μ m 7.6cm x 7.6cm Square	Mylar 6.0 μ m Cut Sheets 1Set=500 Pieces
TF-160-355	Thin Film Mylar 6.0 μ m D8.9cm Circle	Mylar 6.0 μ m D8.9cm 1Set=500 Pieces
TF-160-4	Thin Film Mylar 6.0 μ m 10.16cm x 91.4m Roll	Mylar 6.0 μ m 10.16cm x 91.4m Roll

TF-212	Thin Film Polypropylene 12.0µm 7.6cm x 91.4m Roll	Polypropylene 12.0µm 7.6cm x 91.4m Roll
TF-212-345	Thin Film Polypropylene 12.0µm 7.6cm x 7.6cm Square	Polypropylene 12µm Cut Sheets 1Set=500 Pieces
TF-240	Thin Film Polypropylene 4.0µm 7.6cm x 91.4m Roll	Polypropylene 4.0µm 7.6cm x 91.4m Roll
TF-240-255	Thin Film Polypropylene 4.0µm D6.4cm Circle	Polypropylene 4.0µm D6.4cm 1Set=500 Pieces
TF-240-30	Thin Film Polypropylene 4.0µm D7.6cm Circle	Polypropylene 4.0µm D7.6cm 1Set=100 Pieces
TF-240-345	Thin Film Polypropylene 4.0µm 7.6cm x 7.6cm Square	Polypropylene 4µm Cut Sheets 1Set=500 Pieces
TF-240-355	Thin Film Polypropylene 4.0µm D8.9cm Circle	Polypropylene 4.0µm D8.9cm 1Set=500 Pieces
TF-240-60	Thin Film Polypropylene 4.0µm D6.4cm Circle	Polypropylene 4.0µm D6.4cm 1Set=60 Pieces
TF-250	Thin Film Polypropylene 5.0µm 7.6cm x 91.4m Roll	Polypropylene 5.0µm 7.6cm x 91.4m Roll
TF-250-345	Thin Film Polypropylene 5.0µm 7.6cm x 7.6cm Square	Polypropylene 5.0µm Cut Sheets 1Set=500 Pieces
TF-260	Thin Film Polypropylene 6.0µm 7.6cm x 91.4m Roll	Polypropylene 6.0µm 7.6cm x 91.4m Roll
TF-260-255	Thin Film Polypropylene 6.0µm D6.4cm Circle	Polypropylene 6.0µm D6.4cm 1Set=500 Pieces
TF-260-345	Thin Film Polypropylene 6.0µm 7.6cm x 7.6cm Square	Polypropylene 6.0µm Cut Sheets 1Set=500 Pieces
TF-260-355	Thin Film Polypropylene 6.0µm D8.9cm Circle	Polypropylene 6.0µm D8.9cm 1Set=500 Pieces
TF-475	Thin Film Kapton 7.5µm 7.6cm x 15.2m Roll	Kapton 7.5µm 7.6cm x 15.2m Roll
TF-475-345	Thin Film Kapton 7.5µm 7.6cm x 7.6cm Square	Kapton 7.5µm Cut Sheets 1Set=500 Pieces
TF-500	Thin Film XRF Tape 25mm x 66m Roll	XRF Tape 25mm x 66m Roll
TF-LS-240-255	Thin Film Polypropylene Low Sulfur 4.0µm D6.4cm Circle	Polypropylene Low Sulfur 4.0µm D6.4cm 1Set=500 Pieces

Couppelles Aluminium pour pastillage

AC-530-1	Alu cups straight wall 32mm	Aluminium cups straight wall for pellets 32mm set=1000pcs
AC-532-1	Alu cups taper wall 32mm	Aluminium cups taper wall for pellets 32mm set=1000pcs
AC-540-6	Alu cups taper wall 40mm	Aluminium cups taper wall for pellets 40mm set=600pcs
AC-550-6	Alu cups straight wall 40mm	Aluminium cups straight wall for pellets 40mm set=600pcs

Liants pour pastillage

123	Boreox binder for xrf tablet 0.33g	Boreox binder for xrf tablet 0.33g, diameter 10mm, thickness approx. 5,5mm
1234	Boreox binder for xrf tablet 0.33g	Boreox binder for xrf tablet 0.33g, diameter 10mm, thickness approx. 5,5mm
AC-538-6	Alu cups straight wall 40mm	Aluminium cups straight wall for pellets 40mm set=600pcs
BM-0002-1	Cereox 250g	CEREOX Licowax C Micropowder (old name: Hoechstwax) Wax binder for xrf 250g powder
BM-0002-1B	Cereox 20kg	CEREOX Licowax C Micropowder (old name: Hoechstwax) Wax binder for xrf 20kg powder
BM-0003-1	Cellulose tablet 1kg	Cellulose C20 grinding aid tablets a 0,33g, 1kg
BM-0003-20	Celleox 20kg	Celleox binder for XRF on cellulose basis 20kg
BM-0003-250	Celleox 250g	Celleox binder for XRF on cellulose basis 250g
BM-0003-5	Cellulose tablet 5kg	Cellulose C20 grinding aid tablets a 0,33g, 5kg
BM-0008	Boreox 250g	Boreox, binder for XRF 250g (replacement for boric acid)
BM-0008-1K	Boreox binder for xrf tablet 0.33g	Boreox binder for xrf tablet 0.33g, diameter 10mm, thickness approx. 5,5mm
BM-0008-20K	Boreox binder for xrf tablet 0.33g	Boreox binder for xrf tablet 0.33g, diameter 10mm, thickness approx. 5,5mm
BM-0008-5	Boreox 5kg	Boreox, binder for XRF 5kg (replacement for boric acid)
BM-0008-5K	Boreox binder for xrf tablet 0.33g	Boreox binder for xrf tablet 0.33g, diameter 10mm, thickness approx. 5,5mm

BM-0008-B	Boreox 20kg	BOREOX, binder for XRF, Quantity 20kg (replacement for boric acid)
BM-0008-pellet_34mm	Boreox pellet	Prepared Boreox pellet with 1mm deepening from filling with fat or cream for fatanalysis according to DIN51829, Diameter 34mm, 100 pcs.
BM-0008-pellet_40mm	Boreox pellet	Boreox pellet with 1mm deepening from filling with fat or cream for fatanalysis according to DIN51829, Diameter 40mm,100 pcs.
BM-PXR-225	Multimix 125g	Multimix grinding aid tablets a 0,25g, 125g
BM-PXR-225-1K	Multimix 1kg	Multimix grinding aid tablets a 0,25g, 1kg
BM-PXR-225-25K-Sale	Multimix 25kg	Multimix grinding aid tablets a 0,25g, 25kg
BM-PXR-250	Multimix 250g	Multimix grinding aid tablets a 0,5g, 250g
MU-K-Mixer_50Hz	Mixer to mix powder with wax (220V)	Mixer 220 V 50Hz to mix powder with wax (does not include holding clip for mixing beakers) color: white
MU-K-Mixer_60Hz	Mixer to mix powder with wax (110V)	Mixer 110 V 60Hz to mix powder with wax (does not include holding clip for mixing beakers) color: blue
MU-RMM400	Mixer Mill MM 400 for 100-240V, 50/60Hz	Mixer Mill MM 400 for 100-240V, 50/60Hz
MU-XRF-Mixing_Set_50Hz	XRF-Mixing-Set (220V)	XRF Mixing Starter Kit (Mixer, Holding Clip, non-sliding-mat, 1 tweezer, 5 O-rings, 100 containers, 100 balls, 1 bottle wax)
MU-XRF-Mixing-Set_60Hz	XRF-Mixing-Set (110V)	XRF Mixing Starter Kit (Mixer, Holding Clip, non-sliding-mat, 1 tweezer, 5 O-rings, 100 containers, 100 balls, 1 bottle wax)
TF-PF-32500	Pellet film 32mm	Pellet film 32mm 1 set = 500pcs
TF-PF-40500	Pellet film 40mm	Pellet film 40mm 1 set = 500pcs
WZ-0004	Handling tool for glass beads and pellets	Handling tool for glass beads and pellets

Presses manuel et automatique

PR-25-HD	Vaneox 25t manual	Vaneox -25t manual press
PR-25A-HD	Vaneox 25t automatic	Vaneox 25t automatic press
PR-40-HD	Vaneox 40t automatic	Vaneox 40t automatic press
PR-CD-6032	Die 32mm	Evacuable die complete 32mm

Moules de compression

PR-CD-BA32	Filling tool 32mm	Filling tool for BOREOX (former boric acid) used with a 32mm die
PR-CD-6040	Die 40mm	Evacuable die complete 40mm
PR-CD-6008-6020	Die 8 to 20mm	Evacuable die complete. Diameter from 8 to 20mm according customer specification.
PR-CD-6021-6031	Die 21 to 31mm	Evacuable die complete. Diameter from 21 to 31mm according customer specification.
PR-CD-6033-6039	Die 33 to 39mm	Evacuable die complete. Diameter from 33 to 39mm according customer specification.
PR-CD-6041-6051	Die 41 to 51mm	Evacuable die complete. Diameter from 41 to 51mm according customer specification.
PR-L33000	Die 13mm	Evacuable die complete 13mm
PR-40-CD32	Die 32mm	Evacuable die complete 32mm
PR-40-CD40	Die 40mm	Evacuable die complete 40mm
PR-40-CD35	Die 35mm	Evacuable die complete 35mm
PR-40-CD35R	Die 35mm for ring	Evacuable die complete 35mm for ring OD 51,5mm ID 35mm
PR-40-CD32R	Die 32mm for ring	Evacuable die complete 32mm for ring OD 40mm ID 32mm
PR-CD-6032R	Die 32mm for ring	Evacuable die complete 32mm for ring OD 40mm ID 32mm
PR-CD-6035R_a	Die 35mm for ring	Evacuable die complete 35mm for ring OD 51,5mm ID 35mm

PR-CD-6035R_b	Die 35mm for ring	Evacuable die complete 35mm for ring OD 40mm ID 35mm
PR-CD-BA40	Filling tool 40mm	Filling tool for BOREOX (former boric acid) used with a 40mm die

Machines de fusion gaz

VU-0002-2MA	Vulcan 2MA	Vulcan 2 stations microprocessor + automatic regulation
VU-table	Table for Vulcan Fusion machine and Vaneox Presses	Special table with integrated drawer for Vulcan Fusion machine and Vaneox presses
VU-AB0002-ab	Hood for Vulcan (built-on)	Hood (built-on) for VULCAN without fan Please ask for installation requirements.
VU-0004-6MA	Vulcan 6MA	Vulcan 6 stations microprocessor + automatic regulation
VU-0004-4MA	Vulcan 4MA	Vulcan 4 stations microprocessor + automatic regulation
VU-AB0002-zr	Hood for Vulcan (fragmented)	Hood (fragmented) for VULCAN without fan Please ask for installation requirements.
VU-AB0003-50Hz	High duty radial ventilator 220V	High duty radial ventilator.220V/50Hz suitable for hood VU-AB0002 for fusion machine VULCAN
VU-AB0003-60Hz	High duty radial ventilator 110V	High duty radial ventilator 110V/60Hz suitable for hood VU-AB0002 for fusion machine VULCAN
VU-0006-22A	Vulcan 22MA	Automatic Fusionmachine VULCAN with 2 stations XRF and 2 stations ICP microprocessor control.
VU-0006-21A	Vulcan 21MA	Automatic Fusionmachine VULCAN with 2 stations XRF and 1 station ICP microprocessor control.
VU-0010	Mould feed for VULCAN fusion machine	Mould feed as additional option for VULCAN fusion machine type MA
VU-0011	Mould feed and control for VULCAN fusion machine	Mould feed and control as additional option for VULCAN fusion machine type MA
VU-0016	Zircon crucible for peroxide fusion	Zircon crucible for peroxide fusion
VU-0019	Adapter for VULCAN for zirconia crucibles for peroxide fusion	Adapter for VULCAN fusion machine for 2 - 6 zirconia crucibles for peroxide fusion

VU-VU	Vulcan 22MA	Vulcan 4 Stationen Mikroprozessor + automatische Regelung
VU-0006-6MA-ICP	Vulcan with 6 stations ICP	Automatic Fusionmaschine VULCAN with 6 stations ICP, mikroprozessor control and automatic temperature adjustment
	Vulcan 4MA	Vulcan 4 stations mikroprozessor + automatic regulation
VU-0006-Peroxid	Vulcan 6MA for peroxid fusion	Vulcan fusion system with 6 stations for peroxide fusion
VU-0013	Holder for crucible lid for VULCAN	Holder for crucible lid for VULCAN fusion machine type MA (lids for crucibles have to be purchased separately)

Machines de fusion électrique

VI-0001-1	fully automatic electrical fusion machine for 1 sample, sample holder, protective housing	fully automatic electrical fusion machine for 1 sample, sample holder, protective housing
VI-0001-4	fully automatic electrical fusion machine for 4 samples, incl. autosampler and 4 sample holders, protective housing	fully automatic electrical fusion machine for 4 samples, incl. autosampler and 4 sample holders, protective housing
VI-0001-6	fully automatic electrical fusion machine for 6 samples, incl. autosampler and 6 sample holders, protective housing	fully automatic electrical fusion machine for 6 samples, incl. autosampler and 6 sample holders, protective housing
VI-0001-2	fully automatic electrical fusion machine for 2 samples, incl. autosampler and 2 sample holders, protective housing	fully automatic electrical fusion machine for 2 samples, incl. autosampler and 2 sample holders, protective housing
VI-0002	Compressor for electrical fusion machine	Compressor for compressed air for cooling with automatic dehydration

VI-0003	Exhaust hood for electrical fusion machine	Exhaust hood (only flue) for electrical fusion machine with connection for aspiration
VI-0004	Laboratory bench for electrical fusion machine	Laboratory bench for safety hood for electrical fusion machine
,	fully automatic electrical fusion machine for 2 samples	electrical fusion machine 2 stations microprocessor + automatic regulation
VI-0005	Installation and training for electrical fusion machine	Installation and training for electrical fully automatic electrical fusion machine, 2 days
1	LABORATORY BENCH for electrical fusion machine	Laboratory bench for safety hood for electrical fusion machine
VI-0006	Extension module for handling from 2 up to 4 samples	Extension module for automatic handling of 2 to up 4 samples incl. 2 sample holders
VI-0007	Extension module for handling from 2 up to 6 samples	Extension module for automatic handling of 2 to up to 6 samples incl. 4 sample holders
VI-0100	Sample holder bottom for electrical fusion machine	Sample holder bottom for electrical fusion machine
VI-0101	Adapter for moulds for electrical fusion machine	Adapter for moulds 29/31 mm and 32/34 mm for electrical fusion machine
VI-0102	Sample holder, ceramic crucible holder, 1 piece for electrical fusion machine	Sample holder, holder ceramic crucible holder, 1 piece for electrical fusion machine
VI-0103	Furnace ceramic bottom plate for electrical fusion machine	Furnace ceramic bottom plate for electrical fusion machine
VI-0104	Furnace reservoir sample holder of furnace for electrical fusion machine	Furnace reservoir sample holder of furnace for electrical fusion machine
VI-0105	Furnace inner tube for electrical fusion machine	Furnace inner tube for electrical fusion machine

VI-0106	Heating element for electrical fusion machine	Heating element for electrical fusion machine
VI-0107	Supporting bar 8x8x65 mm for electrical fusion machine	Supporting bar 8x8x65 mm for electrical fusion machine
VI-0108	Cover plate 69x67 mm, 5x2 mm for electrical fusion machine	Cover plate 69x67 mm, 5x2 mm for electrical fusion machine
VI-0109	Temperature switch 36TXE21 for electrical fusion machine	Temperature switch 36TXE21 for electrical fusion machine
VI-0110	Sample holder, complete assembled for electrical fusion machine	Sample holder, complete assembled (VI-0100 + VI-0101) for electrical fusion machine
VI-0008	Rear cover for electrical fusion machine	Rear cover, in case that the electrical fusion machine should be positioned not directly at a wall
VI-0009	High duty radial ventilator for electrical fusion machine	High duty radial ventilator for electrical fusion machine 220V/50Hz
VI-0010	High duty radial ventilator for electrical fusion machine	High duty radial ventilator for electrical fusion machine 110V/60Hz
VI-0012	ICP option including stir plate and teflon beaker for 1 sample	ICP option including stir plate and teflon beaker for 1 sample
WZ-0004	Handling tool for glass beads and pellets	Handling tool for glass beads and pellets
MU-XRF-Mixing_Set_50 Hz	XRF-Mixing-Set (220V)	XRF Mixing Starter Kit (Mixer, Holding Clip, non-sliding-mat, 1 tweezer, 5 O-rings, 100 containers, 100 balls, 1 bottle wax)
MU-XRF-Mixing-Set_60Hz	XRF-Mixing-Set (110V)	XRF Mixing Starter Kit (Mixer, Holding Clip, non-sliding-mat, 1 tweezer, 5 O-rings, 100 containers, 100 balls, 1 bottle wax)

Fondants pour la fusion des perles

FS-OT866	Crucible Platinum/Gold for VULCAN	XRF crucible Platinum/Gold 95/5 for Vulcan
FX-X00-2	Lithiummetaborate 99.98%	Di-Lithiummetaborate 99.98% Pack.=2kg
FX-X20-2	Lithiumtetraborate : Lithiummetaborate 20:80 99,98%	Di-Lithiumtetraborate : Lithiummetaborate 20:80 99,98% Pack.=2kg
FX-X35-2	Lithiumtetraborate : Lithiummetaborate 35:3:64.7 99,98%	Di-Lithiumtetraborate : Lithiummetaborate 35:3:64,7 99,98% Pack.=2kg
FX-X50-2	Lithiumtetraborate : Lithiummetaborate 50:50 99,98%	Di-Lithiumtetraborate : Lithiummetaborate 50:50 99,98% Pack.=2kg
FX-X65-2	Lithiumtetraborate : Lithiummetaborate 66:34 99,98%	Di-Lithiumtetraborate : Lithiummetaborate 66:34 99,98% Pack.=2kg
FX-X80-2	Lithiumtetraborate : Lithiummetaborate 80:20 99,98%	Di-Lithiumtetraborate : Lithiummetaborate 80:20 99,98% Pack.=2kg
FX-X100-2	Lithiumtetraborate 99.98%	Di-Lithiumtetraborate 99.98% Pack.=2kg
FX-X80-2-LiBr	Lithiumtetraborate : Lithiummetaborate 80:20 + 0,4% LiBr 99,98%	Di-Lithiumtetraborate : Lithiummetaborate 80:20 + 0,4% LiBr 99,98% Pack.=2kg
FX-X100-25N	Lithiumtetraborate 99.98%	Di-Lithiumtetraborate 99.98% Pack.=25kg drum
FX-X20-25N	Lithiumtetraborate : Lithiummetaborate 20:80 99,98%	Di-Lithiumtetraborate : Lithiummetaborate 20:80 99,98% Pack.=25kg drum
FX-X50-25N	Lithiumtetraborate : Lithiummetaborate 50:50 99,98%	Di-Lithiumtetraborate : Lithiummetaborate 50:50 99,98% Pack.=25kg drum
FX-X65-25N	Lithiumtetraborate : Lithiummetaborate 66:24 99,98%	Di-Lithiumtetraborate : Lithiummetaborate 66:34 99,98% Pack.=25kg drum
FX-ADD3	Single Additive Non Wetting Agence	Single Additive Non Wetting Agence
FX-CAT8	Single catalyst for ferroalloys FeTi, FeCr, FeV	Single catalyst for ferroalloys FeTi, FeCr, FeV
FX-OXY6	Single oxidizer for ferroalloys and metals	Single oxidizer for ferroalloys and metals

FX-OXY7	Single oxidizer for slags, ferroalloys, oxides	Single oxidizer for slags, ferroalloys, oxides
FX-SUL1	Borate flux for sulfide ores	Borate flux for sulfide ores
FX-SUL2	Borate flux for sulfide ores	Borate flux for sulfide ores
FX-SUL3	Internal standard tungsten oxide	Internal standard tungsten oxide
FX-SUL4	Single oxidizer for sulfides	Single oxidizer for sulfides

Verres de calibration

CS-0002-Set	Application package Ash	Application package Ash 'ready to go'
CS-0001-Set	Application package Cement	Application package cement 'ready to go'
CS-0001/7-Set	Application package Cement + Raw Materials	Application package cement+raw materials 'ready to go'
CS-0007-Set	Application package Raw Materials	Application package raw materials 'ready to go'
CS-0008-Set	Application package Continuous Casting Powder	Application package continous casting powder 'ready to go'
CS-0005-Set	Application package Raw Professional	Application package raw professional 'ready to go'
CS-0023-Set	Application package Slags+Filter Dust	Application package Slags+Filter Dust 'ready to go'
CS-0022-Set	Application package raw material for glass industry	Application package Raw Materials for Glass 'ready to go'
CS-0016-Set	Application package glass	Application package Glass 'ready to go'
CS-0017-Set	Application package Slags	Application package Slags 'ready to go'