

Meinhard

Nébuliseurs et Kits d'introductions inertes

Choisissez votre Gamme

[PFA Micro-Flow Nebulizer](#)

[PFA-ST Nebulizers](#)

[PolyPro-ST Nebulizers](#)

Vous serez automatiquement redirigé à la section souhaitée

Thermo iCap 6000 Series

Thermo	Analab	ESI	Burgener	Meinhard	Description
Glass (8423120514 31)	TR-30-K2			TR-30-K2	Type K glass nebulizer, lapped flush nozzle, high solids, recessed (0.5 mm) capillary, nominal uptake 2mL/min, 0.7 mL/min @30psi.
Glass	TR-50-A2			TR-50-A2	Type A glass nebulizer, lapped flush nozzle, low cost, high sensitivity, nominal uptake 2mL/min, 1L/min @ 50psi
Glass (8423120513 31)	TR-50-C2			TR-50-C2	Type C glass nebulizer, polished, recessed (0.5mm) capillary, high solids, nominal uptake 2mL/min, 1L/min @ 50psi
Glass (8423120513 21)	MP080963			MP080963	V-Groove nebulizer, Noordermeer, 6mm od, 30 psi
PFA	PFA-ST Micro-flow	PFA-ST Micro-flow		ME2040	MicroFlow PFA-ST self-aspirating nebulizer with external ¼-28 threaded connector. Includes one 400µl/min capillary.
	(inclus 1 capillaire au choix)	(inclus 1 capillaire au choix)			
Poly-Pro	ES-4040 Polypro ST MF	ES-4040 Polypro ST MF		ME4040	PolyPro-ST self-aspirating nebulizer with external ¼-28 threaded connector. Includes one 400µl/min capillary.
	(inclus 1 capillaire)				
PEEK (8423120514 71)	PMM4000		PMM4000	PMM4000	PEEK MiraMist high solids nebulizer. Includes an argon line and sample capillary.
PEEK	AM5000		AM5000	AM5000	PEEK AriMist low flow medium solids. Includes an argon line and sample capillary.
Cleaner	PRL-02			PRL-02	Cleaning tool for all quartz and glass concentric nebulizers
Cleaner	ME2046- 0001			ME2046- 0001	Nebulizer Obstruction Removal Kit for ST Nebulizers

Thermo XSeries 2 ICs

Thermo	Analab			Meinhard	Description
Quartz	SILQ-50-A0.5			SILQ-50-A0.5	Plus Series ultra high purity quartz concentric nebulizer, Type A. 1 L/min and 0.5 mL/min @ 50 psig (3.4 bar), 0.7 mm (0.030") inlet for low internal volume.
Quartz	TQP-50-A0.5			TQP-50-A0.5	Plus Series quartz concentric nebulizer, Type A. 1 L/min and 0.5 mL/min @ 50 psig (3.4 bar). 0.7 mm (0.030") inlet for low internal volume.
Quartz	TQ-30-A0.5			TQ-30-A0.5	Type A concentric quartz nebulizer, lapped fl ush nozzle, low background, high sensitivity, nominal uptake 0.5 mL/min, 1 L/min @ 30 psi
Quartz	TQ-30-A1			TQ-30-A1	Type A quartz nebulizer, lapped fl ush nozzle, low background, high sensitivity, nominal uptake 1 mL/min, 1 L/min @ 30 psi (1110720)
Quartz	TQ-30-A3			TQ-30-A3	Type A quartz nebulizer, lapped flush nozzle, low background, high sensitivity, nominal uptake 3 mL/min, 1 L/min @ 30 psi
Quartz	TQ-50-A0.5			TQ-50-A0.5	Type A concentric quartz nebulizer, lapped fl ush nozzle, low background, high sensitivity, nominal uptake 0.5 mL/min, 1 L/min @ 50 psi
Glass	TRP-50-A0.5			TRP-50-A0.5	TR+ glass nebulizer, high sensitivity, low internal volume inlet, nominal uptake 0.5 mL/min, 1 L/min @ 50 psi
Glass	TR-30-A0.5			TR-30-A0.5	Type A glass nebulizer, lapped fl ush nozzle, general purpose, nominal uptake 0.5mL/min, 1L/min @ 30 psi
Glass (1201318)	TR-30-A1			TR-30-A1	Type A glass nebulizer, lapped fl ush nozzle, general purpose, nominal uptake 1mL/min, 1L/min @ 30 psi (1201318)
Glass (1040894)	TR-30-A3			TR-30-A3	Type A glass nebulizer, lapped fl ush nozzle, low cost, high sensitivity, nominal uptake 3 mL/min, 1 L/min @ 30 psi (1040894)
PFA (1301660)	PFA-ST Micro-flow	PFA-ST Micro-flow		ME2040	MicroFlow PFA-ST self-aspirating nebulizer with external ¼-28 threaded connector. Includes one 400µl/min capillary. (1301660)
	(inclus 1 capillaire au choix)	(inclus 1 capillaire au choix)			
Poly-Pro	ES-4040 Polypro ST MF	ES-4040 Polypro ST MF		ME4040	PolyPro-ST self-aspirating nebulizer with external ¼-28 threaded connector. Includes one 400µl/min capillary.
	(inclus 1 capillaire)	(inclus 1 capillaire)			
PEEK (4600356)	PMM4000		PMM3999	PMM4000	PEEK MiraMist high solids nebulizer. Includes an argon line and sample capillary.

PEEK (4600495)	AM5000		AM4999	AM5000	PEEK AriMist low flow medium solids. Includes an argon line and sample capillary.
Cleaner	PRL-02			PRL-02	Cleaning tool for all quartz and glass concentric nebulizers
Cleaner	ME2046-0001			ME2046-0001	Nebulizer Obstruction Removal Kit for ST Nebulizers
Icap-Q					
1161870	TRP-50-A0.5			TRP-50-A0.5	Plus Series Glass concentric, Type A (1 L/min and 0.5 mL/min @ 50 psig)
1162030	ME2050-52			ME2050-52	PFA-LC with tubing
1317090	ME2040-73			ME2040-73	PFA-ST with tubing
1600165	ME2000			ME2000	PFA Concentric, Low Flow (50 uL/min)
1600342	ME2002			ME2002	PFA Concentric, Low Flow (100 uL/min)
1600427	ME2020			ME2020	PFA Concentric, Low Flow (20 uL/min)
4600356	PPM4000			MBR-Peek MiraMist-T	MiraMist, High Solids (0.4-2.0mL/min) HF Resistant
4600495	AP5000			MBR-Peek AriMist-T	AriMist, Low Flow, High Solids (0.2-2.0mL/min) HF Resistant
	PRL-02			PRL-02	Cleaning tool for all quartz and glass concentric nebulizers
Perkin Nexion					
N8145011	TQ-50-A0.5				Quartz concentric, Type A, 1.0 L/min (0.5mL/min, 50 psig)
N8145012	TR-50-C0.5				Glass concentric, Type C, 1.0 L/min (0.5mL/min, 50 psig)
N8122461	TQP-50-A0.5				Quartz Plus - low dead volume, high sensitivity 0.5mL/min (1 L/min, 50 psig)
N0777375	ME2040-7000				Apex-ST PFA MicroFlow Nebulizer
N0777376	ME2002-7000				High temperature APEX-100 PFA Nebulizer
N8145101	ME-ST3-84				PFA ST3 MicroFlow Nebulizer
	PRL-02				Cleaning tool for all quartz and glass concentric nebulizers
Agilent Vista 720/725 ICP-OES Series					
Quartz	SILQ-50-A0.5				Plus Series ultra high purity quartz concentric nebulizer, Type A. 1 L/min and 0.5 mL/min @ 50 psig (3.4 bar), 0.7 mm (0.030") inlet for low internal volume.
Quartz	TQP-50-A0.5				Plus Series quartz concentric nebulizer, Type A. 1 L/min and 0.5 mL/min @ 50 psig (3.4 bar). 0.7 mm (0.030") inlet for low internal volume.

Quartz	TQ-30-A2				Type A quartz nebulizer, lapped flush nozzle, low background, high sensitivity, nominal uptake 2 mL/min, 1 L/min @ 30 psi
Quartz	TQ-30-A2				Type A quartz nebulizer, lapped flush nozzle, low background, high sensitivity, nominal uptake 2 mL/min, 1 L/min @ 30 psi (G1820-65138)
Glass	TRP-50-A0.5				TR+ glass nebulizer, high sensitivity, low internal volume inlet, nominal uptake 0.5 mL/min, 1 L/min @ 50 psi
Glass	TR-30-C2 (2010076500)				Type C glass nebulizer, polished, recessed (0.5mm) capillary, high solids, nominal uptake 2mL/min, 1L/min @ 30 psi
Glass	TR-30-K2 (2010069500)				Type K glass nebulizer, lapped flush nozzle, high solids, recessed (0.5mm) capillary, nominal uptake 2mL/min, 0.7mL/min @30psi.
Glass	TR-30-K3 (2010081600)				Type K glass nebulizer, lapped flush nozzle, high solids, recessed (0.5mm) capillary, nominal uptake 3mL/min, 0.7mL/min @30psi.
PFA	ME2040				MicroFlow PFA-ST self-aspirating nebulizer with external ¼-28 threaded connector. Includes one 400µl/min capillary.
Poly-Pro	ME4040				PolyPro-ST self-aspirating nebulizer with external ¼-28 threaded connector. Includes one 400µl/min capillary.
218080256	MM80256				PEEK Unifit Adapter, 6mm. Glass concentric nebulizer to 12mm inlet port
Cleaner	PRL-02				Cleaning tool for all quartz and glass concentric nebulizers

Agilent 7700 / 8800 Series

Agilent			Meinhard	Description
Quartz	SILQ-50-A0.5		SILQ-50-A0.5	Plus Series ultra high purity quartz concentric nebulizer, Type A. 1 L/min and 0.5 mL/min @ 50 psig (3.4 bar), 0.7 mm (0.030") inlet for low internal volume.
Quartz	TQP-50-A0.5		TQP-50-A0.5	Plus Series quartz concentric nebulizer, Type A. 1 L/min and 0.5 mL/min @ 50 psig (3.4 bar). 0.7 mm (0.030") inlet for low internal volume.
Quartz	TQ-30-A0.5		TQ-30-A0.5	Type A concentric quartz nebulizer, lapped fl ush nozzle, low background, high sensitivity, nominal uptake 0.5 mL/min, 1 L/min @ 30 psi
Quartz	TQ-30-A1		TQ-30-A1	Type A quartz nebulizer, lapped fl ush nozzle, low background, high sensitivity, nominal uptake 1 mL/min, 1 L/min @ 30 psi
Quartz	TQ-30-A2		TQ-30-A2	Type A quartz nebulizer, lapped fl ush nozzle, low background, high sensitivity, nominal uptake 2 mL/min, 1 L/min @ 30 psi (G1820-65138)
Glass	TRP-50-A0.5		TRP-50-A0.5	TR+ glass nebulizer, high sensitivity, low internal volume inlet, nominal uptake 0.5 mL/min, 1 L/min @ 50 psi
Glass	TR-30-A0.5		TR-30-A0.5	Type A glass nebulizer, lapped fl ush nozzle, general purpose, nominal uptake 0.5mL/min, 1L/min @ 30 psi
Glass	TR-30-A1		TR-30-A1	Type A glass nebulizer, lapped fl ush nozzle, general purpose, nominal uptake 1mL/min, 1L/min @ 30 psi (G1820-65030)
PFA	ME2040		ME2040	MicroFlow PFA-ST self-aspirating nebulizer with external ¼-28 threaded connector. Includes one 400µl/min capillary.
Poly-Pro	ME4040		ME4040	PolyPro-ST self-aspirating nebulizer with external ¼-28 threaded connector. Includes one 400µl/min capillary.
Cleaner	PRL-02		PRL-02	Cleaning tool for all quartz and glass concentric nebulizers