



# Pompes péristaltiques et tuyaux

Choisissez votre Gamme

Pompes péristaltiques

[Pompes PERIMAX](#)

[Pompes SYMAX](#)

Tuyaux

[PVC](#)

[Solvent Resistant](#)

[Viton](#)

[Santoprene](#)

[Silicone](#)

Vous serez automatiquement redirigé à la section souhaitée

## PERIMAX® 12

PERIMAX® 12/1 (1 channel) 08-1000

PERIMAX® 12/2 (2 channel) 08-1001

PERIMAX® 12/4 (4 channel) 08-1003

PERIMAX® 12/6 (6 channel) 08-0104

### Option:

Remote Control Gilson compatible (08-0148)

ON / OFF, high speed, clockwise / counter-clockwise, rpm

## PERIMAX® 16 - NON-PULSING

PERIMAX® 16/1 (1 pulse-free channel) 08-1501

PERIMAX® 16/2 (2 pulse-free channels) 08-1502

PERIMAX® 16/3 (3 pulse-free channels) 08-0503

### Option:

Remote Control Gilson compatible (08-0148)

ON / OFF, high speed, clockwise / counter-clockwise, rpm

### Tubings:

PVC tubes (1/2 dozen)






pre-assembled by Y-junction

# SYMAX





No.	Description
30/09/02	Syringe pump Symax incl. adjustable syringe holder and electrical connection kit
10-0103	Syringe pump Symax (additional channel) incl. adjustable syringe holder
<b>Equipment:</b>	
Syringe	
10-0010	Polypropylene syringe, 50 ml
10-0011	Polypropylene syringe, 5ml
10-0013	PTFE syringe luer 10ml
10-0015	Glass syringe luer, 1ml
10-0016	Glass syringe luer, 5ml
10-0017	Glass syringe luer, 10ml
10-0018	Glass syringe luer, 20ml
10-0019	Glass syringe luer, 30ml
10-0020	Glass syringe luer, 50ml
Sets of tubing	
10-0082	PTFE set of tubing
10-0081	PVC set of tubing
10-0080	3-way-valve set
10-0060	Electrical connection kit
10-0050	Foot switch

<b>10-0061</b>	Connecting cable, 1m
<b>10-0062</b>	OEM network cable, 2m
<b>10-0063</b>	OEM trigger cable, 2m

## PVC Pump Tubing Standard



					 <b>2 bridges</b> 12 pieces per package	 <b>2 bridges</b> 12 pieces per package	 <b>3 bridges</b> 12 pieces per package	 <b>3 bridges</b> 12 pieces per package	 <b>3 bridges</b> 12 pieces per package
mm	inch	mm	inch		Part Number	Part Number	Part Number	Part Number	Part Number
0,127	0,0050	2,007	0,079	orange/black	38-0005	381-0005	381-0205	38-0205	38-0405
0,191	0,0075	2,032	0,080	orange/red	38-0010	381-0010	381-0210	38-0210	38-0410
0,254	0,0100	2,083	0,082	orange/blue	38-0015	381-0015	381-0215	38-0215	38-0415
0,381	0,0150	2,209	0,087	orange/green	38-0020	381-0020	381-0220	38-0220	38-0420
0,508	0,0200	2,337	0,092	orange/yellow	38-0025	381-0025	381-0225	38-0225	38-0425
0,635	0,0250	2,464	0,097	orange/white	38-0030	381-0030	381-0230	38-0230	38-0430
0,762	0,0300	2,438	0,096	black/black	38-0032	381-0032	381-0232	38-0232	38-0432
0,889	0,0350	2,565	0,101	orange/orange	38-0035	381-0035	381-0235	38-0235	38-0435
1,016	0,0400	2,692	0,106	white/white	38-0040	381-0040	381-0240	38-0240	38-0440
1,143	0,0450	2,819	0,111	red/red	38-0043	381-0043	381-0243	38-0243	38-0443
1,295	0,0510	2,972	0,117	grey/grey	38-0046	381-0046	381-0246	38-0246	38-0446
1,422	0,0560	3,099	0,112	yellow/yellow	38-0050	381-0050	381-0250	38-0250	38-0450
1,524	0,0600	3,200	0,126	yellow/blue	38-0052	381-0052	381-0252	38-0252	38-0452
1,651	0,0650	3,327	0,131	blue/blue	38-0055	381-0055	381-0255	38-0255	38-0455
1,750	0,0690	3,430	0,135	blue/green	38-0058	381-0058	381-0258	38-0258	38-0458
1,854	0,0730	3,531	0,139	green/green	38-0060	381-0060	381-0260	38-0260	38-0460
2,057	0,0810	3,734	0,147	purple/purple	38-0065	381-0065	381-0265	38-0265	38-0465
2,286	0,0900	3,962	0,156	purple/black	38-0067	381-0067	381-0267	38-0267	38-0467
2,540	0,1000	4,216	0,166	purple/orange	38-0070	381-0070	381-0270	38-0270	38-0470
2,794	0,1100	4,470	0,176	purple/white	38-0078	381-0078	381-0278	38-0278	38-0478
3,175	0,1250	4,851	0,191	black/white	38-0080	381-0080	381-0280	38-0280	38-0480

## PVC Pump Tubing Standard with Flared End






					 <b>2 bridges</b> 12 pieces per package	 <b>2 bridges</b> 12 pieces per package	 <b>3 bridges</b> 12 pieces per package	 <b>3 bridges</b> 12 pieces per package
Inside Diameter:		Outside Diameter:		Colour Code:	Part Number	Part Number	Part Number	Part Number
mm	inch	mm	inch		Part Number	Part Number	Part Number	Part Number
0,127	0,0050	2,007	0,079	orange/black	38-0005F	381-0005F	381-0205F	38-0205F
0,191	0,0075	2,032	0,080	orange/red	38-0010F	381-0010F	381-0210F	38-0210F
0,254	0,0100	2,083	0,082	orange/blue	38-0015F	381-0015F	381-0215F	38-0215F
0,381	0,0150	2,209	0,087	orange/green	38-0020F	381-0020F	381-0220F	38-0220F
0,508	0,0200	2,337	0,092	orange/yellow	38-0025F	381-0025F	381-0225F	38-0225F
0,635	0,0250	2,464	0,097	orange/white	38-0030F	381-0030F	381-0230F	38-0230F

0,762	0,0300	2,438	0,096	black/black	38-0032F	381-0032F	381-0232F	38-0232F
0,889	0,0350	2,565	0,101	orange/orange	38-0035F	381-0035F	381-0235F	38-0235F
1,016	0,0400	2,692	0,106	white/white	38-0040F	381-0040F	381-0240F	38-0240F



## PVC Pump Tubing Standard

					by the meter	by the mete
						
Inside Diameter:		Outside Diameter:		Colour Code:	3 m (118.11 inch) per package	10 m (393.70 inch) per package
mm	inch	mm	inch		Part Number	Part Number
0,127	0,0050	2,007	0,079	orange/black	38-0105	38-0505
0,191	0,0075	2,032	0,080	orange/red	38-0110	38-0510
0,254	0,0100	2,083	0,082	orange/blue	38-0115	38-0515
0,381	0,0150	2,209	0,087	orange/green	38-0120	38-0520
0,508	0,0200	2,337	0,092	orange/yellow	38-0125	38-0525
0,635	0,0250	2,464	0,097	orange/white	38-0130	38-0530
0,762	0,0300	2,438	0,096	black/black	38-0132	38-0532
0,889	0,0350	2,565	0,101	orange/orange	38-0135	38-0535
1,016	0,0400	2,692	0,106	white/white	38-0140	38-0540
1,143	0,0450	2,819	0,111	red/red	38-0143	38-0543
1,295	0,0510	2,972	0,117	grey/grey	38-0146	38-0546
1,422	0,0560	3,099	0,112	yellow/yellow	38-0150	38-0550
1,524	0,0600	3,200	0,126	yellow/blue	38-0152	38-0552
1,651	0,0650	3,327	0,131	blue/blue	38-0155	38-0555
1,750	0,0690	3,430	0,135	blue/green	38-0158	38-0558
1,854	0,0730	3,531	0,139	green/green	38-0160	38-0560
2,057	0,0810	3,734	0,147	purple/purple	38-0165	38-0565
2,286	0,0900	3,962	0,156	purple/black	38-0167	38-0567
2,540	0,1000	4,216	0,166	purple/orange	38-0170	38-0570
2,794	0,1100	4,470	0,176	purple/white	38-0178	38-0578
3,175	0,1250	4,851	0,191	black/white	38-0180	38-0580

## PU Longlife Solvent Resistant Flex

					2 bridges 	2 bridges 	3 bridges 	3 bridges 	3 bridges 
Inside Diameter:		Outside Diameter:		Colour Code:	12 pieces per package	12 pieces per package	12 pieces per package	12 pieces per package	12 pieces per package
mm	inch	mm	inch		Part Number	Part Number	Part Number	Part Number	Part Number
0,127	0,0050	2,007	0,079	orange/black	38-2005	381-2005	381-2205	38-2205	38-2405
0,191	0,0075	2,032	0,080	orange/red	38-2010	381-2010	381-2210	38-2210	38-2410
0,254	0,0100	2,083	0,082	orange/blue	38-2015	381-2015	381-2215	38-2215	38-2415
0,381	0,0150	2,209	0,087	orange/green	38-2020	381-2020	381-2220	38-2220	38-2420
0,508	0,0200	2,337	0,092	orange/yellow	38-2025	381-2025	381-2225	38-2225	38-2425
0,635	0,0250	2,464	0,097	orange/white	38-2030	381-2030	381-2230	38-2230	38-2430
0,762	0,0300	2,438	0,096	black/black	38-2032	381-2032	381-2232	38-2232	38-2432
0,889	0,0350	2,565	0,101	orange/orange	38-2035	381-2035	381-2235	38-2235	38-2435
1,016	0,0400	2,692	0,106	white/white	38-2040	381-2040	381-2240	38-2240	38-2440
1,143	0,0450	2,819	0,111	red/red	38-2043	381-2043	381-2243	38-2243	38-2443
1,295	0,0510	2,972	0,117	grey/grey	38-2046	381-2046	381-2246	38-2246	38-2446
1,422	0,0560	3,099	0,112	yellow/yellow	38-2050	381-2050	381-2250	38-2250	38-2450
1,524	0,0600	3,200	0,126	yellow/blue	38-2052	381-2052	381-2252	38-2252	38-2452
1,651	0,0650	3,327	0,131	blue/blue	38-2055	381-2055	381-2255	38-2255	38-2455
1,750	0,0690	3,430	0,135	blue/green	38-2058	381-2058	381-2258	38-2258	38-2458
1,854	0,0730	3,531	0,139	green/green	38-2060	381-2060	381-2260	38-2260	38-2460
2,057	0,0810	3,734	0,147	purple/purple	38-2065	381-2065	381-2265	38-2265	38-2465
2,286	0,0900	3,962	0,156	purple/black	38-2067	381-2067	381-2267	38-2267	38-2467
2,540	0,1000	4,216	0,166	purple/orange	38-2070	381-2070	381-2270	38-2270	38-2470
2,794	0,1100	4,470	0,176	purple/white	38-2078	381-2078	381-2278	38-2278	38-2478
3,175	0,1250	4,851	0,191	black/white	38-2080	381-2080	381-2280	38-2280	38-2480






## PU Longlife Solvent Resistant Flex

					by the meter 	by the meter 
Inside Diameter:		Outside Diameter:		Colour Code:	3 m (118.11 inch) per package	10 m (393.70 inch) per package
mm	inch	mm	inch		Part Number	Part Number
0,127	0,0050	2,007	0,079	orange/black	38-2105	38-2505
0,191	0,0075	2,032	0,080	orange/red	38-2110	38-2510
0,254	0,0100	2,083	0,082	orange/blue	38-2115	38-2515
0,381	0,0150	2,209	0,087	orange/green	38-2120	38-2520
0,508	0,0200	2,337	0,092	orange/yellow	38-2125	38-2525
0,635	0,0250	2,464	0,097	orange/white	38-2130	38-2530
0,762	0,0300	2,438	0,096	black/black	38-2132	38-2532




0,889	0,0350	2,565	0,101	orange/orange	38-2135	38-2535
1,016	0,0400	2,692	0,106	white/white	38-2140	38-2540
1,143	0,0450	2,819	0,111	red/red	38-2143	38-2543
1,295	0,0510	2,972	0,117	grey/grey	38-2146	38-2546
1,422	0,0560	3,099	0,112	yellow/yellow	38-2150	38-2550
1,524	0,0600	3,200	0,126	yellow/blue	38-2152	38-2552
1,651	0,0650	3,327	0,131	blue/blue	38-2155	38-2555
1,750	0,0690	3,430	0,135	blue/green	38-2158	38-2558
1,854	0,0730	3,531	0,139	green/green	38-2160	38-2560
2,057	0,0810	3,734	0,147	purple/purple	38-2165	38-2565
2,286	0,0900	3,962	0,156	purple/black	38-2167	38-2567
2,540	0,1000	4,216	0,166	purple/orange	38-2170	38-2570
2,794	0,1100	4,470	0,176	purple/white	38-2178	38-2578
3,175	0,1250	4,851	0,191	black/white	38-2180	38-2580



## VITON® Pump Tubing Acid Resistant Flex






					2 bridges	2 bridges	3 bridges	3 bridges	3 bridges
									
Inside Diameter:		Outside Diameter:		Colour Code:	12 pieces per package	12 pieces per package	12 pieces per package	12 pieces per package	12 pieces per package
mm	inch	mm	inch		Part Number	Part Number	Part Number	Part Number	Part Number
0,508	0,0200	2,337	0,092	orange/yellow	38-5625	381-5625	381-5225	38-5225	38-5425
0,635	0,0250	2,464	0,097	orange/white	38-5630	381-5630	381-5230	38-5230	38-5430
0,762	0,0300	2,438	0,096	black/black	38-5632	381-5632	381-5232	38-5232	38-5432
0,889	0,0350	2,565	0,101	orange/orange	38-5635	381-5635	381-5235	38-5235	38-5435
1,016	0,0400	2,692	0,106	white/white	38-5640	381-5640	381-5240	38-5240	38-5440
1,143	0,0450	2,819	0,111	red/red	38-5643	381-5643	381-5243	38-5243	38-5443
1,295	0,0510	2,972	0,117	grey/grey	38-5646	381-5646	381-5246	38-5246	38-5446
1,422	0,0560	3,099	0,112	yellow/yellow	38-5650	381-5650	381-5250	38-5250	38-5450
1,524	0,0600	3,200	0,126	yellow/blue	38-5652	381-5652	381-5252	38-5252	38-5452
1,651	0,0650	3,327	0,131	blue/blue	38-5655	381-5655	381-5255	38-5255	38-5455
1,750	0,0690	3,430	0,135	blue/green	38-5658	381-5658	381-5258	38-5258	38-5458
1,854	0,0730	3,531	0,139	green/green	38-5660	381-5660	381-5260	38-5260	38-5460
2,057	0,0810	3,734	0,147	purple/purple	38-5665	381-5665	381-5265	38-5265	38-5465
2,286	0,0900	3,962	0,156	purple/black	38-5667	381-5667	381-5267	38-5267	38-5467
2,540	0,1000	4,216	0,166	purple/orange	38-5670	381-5670	381-5270	38-5270	38-5470
2,794	0,1100	4,470	0,176	purple/white	38-5678	381-5678	381-5278	38-5278	38-5478

## VITON® Pump Tubing Acid Resistant Flex



					2 bridges (short)	by the meter	by the meter
							
Inside Diameter:		Outside Diameter:		Colour Code:	12 pieces per package	3 m (118.11 inch) per package	10 m (393.70 inch) per package
mm	inch	mm	inch		Part Number	Part Number	Part Number
0,508	0,0200	2,337	0,092	orange/yellow	38-5025	38-5125	38-5525
0,635	0,0250	2,464	0,097	orange/white	38-5030	38-5130	38-5530
0,762	0,0300	2,438	0,096	black/black	38-5032	38-5132	38-5532
0,889	0,0350	2,565	0,101	orange/orange	38-5035	38-5135	38-5535
1,016	0,0400	2,692	0,106	white/white	38-5040	38-5140	38-5540
1,143	0,0450	2,819	0,111	red/red	38-5043	38-5143	38-5543
1,295	0,0510	2,972	0,117	grey/grey	38-5046	38-5146	38-5546
1,422	0,0560	3,099	0,112	yellow/yellow	38-5050	38-5150	38-5550
1,524	0,0600	3,200	0,126	yellow/blue	38-5052	38-5152	38-5552
1,651	0,0650	3,327	0,131	blue/blue	38-5055	38-5155	38-5555
1,750	0,0690	3,430	0,135	blue/green	38-5058	38-5158	38-5558

1,854	0,0730	3,531	0,139	green/green	38-5060	38-5160	38-5560
2,057	0,0810	3,734	0,147	purple/purple	38-5065	38-5165	38-5565
2,286	0,0900	3,962	0,156	purple/black	38-5067	38-5167	38-5567
2,540	0,1000	4,216	0,166	purple/orange	38-5070	38-5170	38-5570
2,794	0,1100	4,470	0,176	purple/white	38-5078	38-5178	38-5578

**Santoprene® TPE Pump Tubing comparable to Pharmed® / Mediprene®**






					2 bridges 	2 bridges 	3 bridges 	3 bridges 	3 bridges 
Inside Diameter:		Outside Diameter:		Colour Code:	12 pieces per package	12 pieces per package	12 pieces per package	12 pieces per package	12 pieces per package
mm	inch	mm	inch		Part Number	Part Number	Part Number	Part Number	Part Number
0,254	0,0100	2,083	0,082	orange/blue	38-6015	381-6015	381-6215	38-6215	38-6415
0,381	0,0150	2,209	0,087	orange/green	38-6020	381-6020	381-6220	38-6220	38-6420
0,508	0,0200	2,337	0,092	orange/yellow	38-6025	381-6025	381-6225	38-6225	38-6425
0,635	0,0250	2,464	0,097	orange/white	38-6030	381-6030	381-6230	38-6230	38-6430
0,762	0,0300	2,438	0,096	black/black	38-6032	381-6032	381-6232	38-6232	38-6432
0,889	0,0350	2,565	0,101	orange/orange	38-6035	381-6035	381-6235	38-6235	38-6435
1,016	0,0400	2,692	0,106	white/white	38-6040	381-6040	381-6240	38-6240	38-6440
1,143	0,0450	2,819	0,111	red/red	38-6043	381-6043	381-6243	38-6243	38-6443
1,295	0,0510	2,972	0,117	grey/grey	38-6046	381-6046	381-6246	38-6246	38-6446
1,422	0,0560	3,099	0,112	yellow/yellow	38-6050	381-6050	381-6250	38-6250	38-6450
1,524	0,0600	3,200	0,126	yellow/blue	38-6052	381-6052	381-6252	38-6252	38-6452
1,651	0,0650	3,327	0,131	blue/blue	38-6055	381-6055	381-6255	38-6255	38-6455
1,750	0,0690	3,430	0,135	blue/green	38-6058	381-6058	381-6258	38-6258	38-6458
1,854	0,0730	3,531	0,139	green/green	38-6060	381-6060	381-6260	38-6260	38-6460
2,057	0,0810	3,734	0,147	purple/purple	38-6065	381-6065	381-6265	38-6265	38-6465
2,286	0,0900	3,962	0,156	purple/black	38-6067	381-6067	381-6267	38-6267	38-6467
2,540	0,1000	4,216	0,166	purple/orange	38-6070	381-6070	381-6270	38-6270	38-6470
2,794	0,1100	4,470	0,176	purple/white	38-6078	381-6078	381-6278	38-6278	38-6478
3,175	0,1250	4,851	0,191	black/white	38-6080	381-6080	381-6280	38-6280	38-6480

**Santoprene® TPE Pump Tubing comparable to Pharmed® / Mediprene®**



					by the meter 	by the meter 
Inside Diameter:		Outside Diameter:		Colour Code:	3 m (118.11 inch) per package	10 m (393.70 inch) per package
mm	inch	mm	inch		Part Number	Part Number
0,254	0,0100	2,083	0,082	orange/blue	38-6115	38-6515
0,381	0,0150	2,209	0,087	orange/green	38-6120	38-6520
0,508	0,0200	2,337	0,092	orange/yellow	38-6125	38-6525
0,635	0,0250	2,464	0,097	orange/white	38-6130	38-6530
0,762	0,0300	2,438	0,096	black/black	38-6132	38-6532
0,889	0,0350	2,565	0,101	orange/orange	38-6135	38-6535
1,016	0,0400	2,692	0,106	white/white	38-6140	38-6540
1,143	0,0450	2,819	0,111	red/red	38-6143	38-6543

1,295	0,0510	2,972	0,117	grey/grey	38-6146	38-6546
1,422	0,0560	3,099	0,112	yellow/yellow	38-6150	38-6550
1,524	0,0600	3,200	0,126	yellow/blue	38-6152	38-6552
1,651	0,0650	3,327	0,131	blue/blue	38-6155	38-6555
1,750	0,0690	3,430	0,135	blue/green	38-6158	38-6558
1,854	0,0730	3,531	0,139	green/green	38-6160	38-6560
2,057	0,0810	3,734	0,147	purple/purple	38-6165	38-6565
2,286	0,0900	3,962	0,156	purple/black	38-6167	38-6567
2,540	0,1000	4,216	0,166	purple/orange	38-6170	38-6570
2,794	0,1100	4,470	0,176	purple/white	38-6178	38-6578
3,175	0,1250	4,851	0,191	black/white	38-6180	38-6580

## Silicone Pump Tubing

					2 bridges	2 bridges	3 bridges	3 bridges	3 bridges
									
Inside Diameter:		Outside Diameter:		Colour Code:	6 pieces per package	6 pieces per package	6 pieces per package	6 pieces per package	6 pieces per package
mm	inch	mm	inch		Part Number	Part Number	Part Number	Part Number	Part Number
0,508	0,0200	2,337	0,092	orange/yellow	38-3025	381-3025	381-3225	38-3225	38-3425
0,635	0,0250	2,464	0,097	orange/white	38-3030	381-3030	381-3230	38-3230	38-3430
0,762	0,0300	2,438	0,096	black/black	38-3032	381-3032	381-3232	38-3232	38-3432
0,889	0,0350	2,565	0,101	orange/orange	38-3035	381-3035	381-3235	38-3235	38-3435
1,016	0,0400	2,692	0,106	white/white	38-3040	381-3040	381-3240	38-3240	38-3440
1,143	0,0450	2,819	0,111	red/red	38-3043	381-3043	381-3243	38-3243	38-3443
1,295	0,0510	2,972	0,117	grey/grey	38-3046	381-3046	381-3246	38-3246	38-3446
1,422	0,0560	3,099	0,112	yellow/yellow	38-3050	381-3050	381-3250	38-3250	38-3450
1,524	0,0600	3,200	0,126	yellow/blue	38-3052	381-3052	381-3252	38-3252	38-3452
1,651	0,0650	3,327	0,131	blue/blue	38-3055	381-3055	381-3255	38-3255	38-3455
1,750	0,0690	3,430	0,135	blue/green	38-3058	381-3058	381-3258	38-3258	38-3458
1,854	0,0730	3,531	0,139	green/green	38-3060	381-3060	381-3260	38-3260	38-3460
2,057	0,0810	3,734	0,147	purple/purple	38-3065	381-3065	381-3265	38-3265	38-3465
2,286	0,0900	3,962	0,156	purple/black	38-3067	381-3067	381-3267	38-3267	38-3467
2,540	0,1000	4,216	0,166	purple/orange	38-3070	381-3070	381-3270	38-3270	38-3470
2,794	0,1100	4,470	0,176	purple/white	38-3078	381-3078	381-3278	38-3278	38-3478

## Silicone Pump Tubing

					by the meter	by the meter
						
					3 m (118.11 inch) per package	10 m (393.70 inch) per package
Inside Diameter:		Outside Diameter:		Colour Code:	Part Number	Part Number
mm	inch	mm	inch			
0,508	0,0200	2,337	0,092	orange/yellow	38-3125	38-3525
0,635	0,0250	2,464	0,097	orange/white	38-3130	38-3530
0,762	0,0300	2,438	0,096	black/black	38-3132	38-3532
0,889	0,0350	2,565	0,101	orange/orange	38-3135	38-3535
1,016	0,0400	2,692	0,106	white/white	38-3140	38-3540
1,143	0,0450	2,819	0,111	red/red	38-3143	38-3543
1,295	0,0510	2,972	0,117	grey/grey	38-3146	38-3546
1,422	0,0560	3,099	0,112	yellow/yellow	38-3150	38-3550
1,524	0,0600	3,200	0,126	yellow/blue	38-3152	38-3552
1,651	0,0650	3,327	0,131	blue/blue	38-3155	38-3555
1,750	0,0690	3,430	0,135	blue/green	38-3158	38-3558

1,854	0,0730	3,531	0,139	green/green	38-3160	38-3560
2,057	0,0810	3,734	0,147	purple/purple	38-3165	38-3565
2,286	0,0900	3,962	0,156	purple/black	38-3167	38-3567
2,540	0,1000	4,216	0,166	purple/orange	38-3170	38-3570
2,794	0,1100	4,470	0,176	purple/white	38-3178	38-3578