

Magnetic Stirrer

Perfect for homogenization of viscous solutions.



Non contractual pictures

Magnetic stirrer allows to homogenize a solution using magnetic force.

This device is mainly used in the following matters:

- prepare a viscous solution
- accelerate dissolution of a solid composite in a solvent
- proceed to a recrystallization

The magnetic stirrer can be used with A4 and A3 hotplates. The screen of the box indicates rotation speed of the magnets (magnets are not sold with the device) and then allows settling the best speed for mineralization.

Important : Temperature of the hotplate cannot exceed 120°C once it is used with a magnetic stirrer.

