

EvapoClean® - Accessories

ELBOW Stand

Non contractual pictures



The “Elbow Stand” accessory was designed to preserve the integrity of your samples after evaporation or your acids after the purification process.

It can be used in the following cases:

- To safely carry samples between the device and the Sorbonne,
- To let stand the elbows before or after used.

Our models can be used with all our elbows (125mL / 25mL / 7mL).

Cales en graphite et adaptateur PTFE

Non contractual pictures



The PFA-coated graphite wedge and PTFE adapter have been designed for users who wish to purify acids in smaller or larger containers than those supplied with the 25ml, 125ml and 500ml devices.


Evapo Cleaning kit

The Cleaning Kit allows you to clean vessel used with your Evapoclean (25ml, 125ml or 500ml).

Each kit includes:

- 3 PFA tanks (30ml or 90ml)
- 3 PTFE sleeves allowing to associate containers that need to be cleaned to tanks.

Working principle



Evapo Cleaning Kit 25ml

3 CleanClever® n°7
3 Vials 30ml

Evapo Cleaning Kit 125ml

3 CleanClever® n°6
3 Jars 90ml

Evapo Cleaning Kit 500ml

3 CleanClever® n°25
3 Jars 90ml