

NEW MARKED VESSEL KIT

VESSEL KIT

As per as our update of the 1st April, we have made some changes related of our Evapoclean and Cleanacids devices. From now on, the device's price will exclude the vessel. In order for you to buy the vessel along with the device, you have discovered new products called "NEW MARKED VESSEL KIT" on the recent pricelist.

These kits are the vessel you can include in the device. Knowing that one Kit include **1 ELBOW + 2 VIALS + 2 CLOSURES + LASER MARKING**.

For example: 1 EV-6.25ml → 6 KIT VESSEL 25ML



LASER MARKING

In order to follow the analysis circuit, to accurately identify and trace all of your labware, we set up the laser marking. The marking will be standard to all of our kit, but you will be able to ask for a customization. The marking will always follow two sequences for each vial in the kit 1A & 1B. If you have 6 Kits, the sequences will end at 6A & 6B.

Standard laser Marking

