


PFA Jars (OD 50 mm)						
Sample jars OD 50mm with thread (without screw closure), inner bottom conical and rounded, outside flat with recess for tool						
50 ml	49	45	52	Total volume 60 ml	T44-405	
100 ml	49	70	52	Total volume 60 ml	T44-410	
180 ml	49	120	52	Total volume 185 ml	T44-418	
Accessories for PAF Jars						
Closure	61	19	52x4	Screw Cap closure for jars T44-4xxfits to wrench tool T44-492	T44-400	
Closure	61	19 + 20	52x4	Screw Cap with Port closure for jars T44-4xx with port tube OD 6 x 20mm	T44-4006	
Wrench Set				Hand Wrench Tool fits in bottom recess of T44- 4xx & rips of closure T44- 400 for tight closing / opening of PFA jars	T44-492	
PFA Bottles						
20 ml	28	60		Wide Neck Bottle with PFA screw capOD 28mm, Height 60mm,Neck ID 16mm	T84-020	
100 ml	48	93		Wide Neck Bottle with PFA screw capOD 48mm, Height 93mm,Neck ID 33mm	T84-100	
PFA Tube						
10 ml	16	109		Sample Tube 16mm without closure Neck ID 13mm. Stoppers and Snap-Caps optional	T43-011	
PFA Marking						



Marks	On request	Laser Marking (serial > 50 pcs.) 4 digits and 2D-codedurable lettering without contamination	T99-XXX	
--------------	------------	---	----------------	---