


ETC EasyTraceCleaner® - Options

OPTION 01 Middle Plate

This middle plate doubles your cleaning capacity by allowing you to wash your vessel above your vials.



Reference	Designation
ETC-O1-L	PTFE MIDDLE PLATE FOR ETC EVO II -L- (Spacer 200mm)
ETC-O1-XL	PTFE MIDDLE PLATE FOR ETC EVO II -XL- (Spacer 250mm)

OPTION 02

Additional Tray

This Tray includes 1 lower plate + 1 handle + 48 Rods.

This option allows you to prepare a second set of Vessel to be washed in advance.




Reference	Designation
ETC-O2-L	ADDITIONNAL PTFE SET: TRAY + HANDLE + 48 RODS (120mm) FOR ETC EVO II L
ETC-O2-XL	ADDITIONNAL PTFE SET: TRAY + HANDLE + 48 RODS (120mm) FOR ETC EVO II XL

OPTION 03

Additional PTFE 48 Rods

Allows you to stabilize your vessel more easily. Transforms a middle plate for closures into a middle plate for vessel





Double rods level only possible with an ETC Size XL

Reference	Designation
ETC-O3-120	ADDITIONNAL PTFE 48 RODS 120mm FOR ETC EVO II
ETC-O3-160	ADDITIONNAL PTFE 48 RODS 160mm FOR ETC EVO II
ETC-O3-210	ADDITIONNAL PTFE 48 RODS 210mm FOR ETC EVO II

OPTION 04/04B

P'TIPS RACK & P'TIPS TRANSFER RACK

Allows you to easily and efficiently clean all your pipette tips in your ETC.




Reference	Designation
ETC-O4-200	P'TIPS RACK 200µl - L115 FOR ETC EVO II
ETC-O4-1000	P'TIPS RACK 1000µl - L115 FOR ETC EVO II
ETC-O4-5000	P'TIPS RACK 5000µl - L115 FOR ETC EVO II
ETC-O4-10000	P'TIPS RACK 10000µl - L115 FOR ETC EVO II
ETC-O4B-T200	P'TIPS TRANSFER RACK 200µL PLATE + RODS 96 POSITIONS FOR ETC EVO II
ETC-O4B-T1000	P'TIPS TRANSFER RACK 1000µL PLATE + RODS 96 POSITIONS FOR ETC EVO II

OPTION 05

Small OD Vessel Holder

This Vessel holder can be used to stabilize your vessel instead of your rods (picture 2). Essential for holding vessel that cannot be held by the PTFE rods (picture 3).

The height of the holder can be adapted according to the size of the vessel to be cleaned.




Reference	Designation
ETC-O5	SMALL OD VESSEL HOLDER FOR ETC EVO II (Spacer 205mm)

OPTION 01/05B

Additional Spacer

Allows you to customize the height of your trays (middle plate or Vessel holder) according to your specific needs.

Contact us to define the height of the Spacer in accordance with your Vessel.



Reference	Designation
ETC-01/05B	CUSTOM SPACER FOR PTFE MIDDLE PLATE OR VESSEL HOLDER FOR ETC EVO II

OPTION 06

Fume Hood for ETC

Helps maintain control of your work environment



Reference	Designation
ETC-06	MINI FUME HOOD FOR ETC EVO II (500X500X600mm)

OPTION 06B

Charcoal Filter

The Charcoal filter system provides optimum safety during your chemical handling



Reference	Designation
ETC-O6B	CHARCOAL FILTER FOR MINI FUME HOOD - ETC EVO II